

TRANSFORMING PATIENT LIVES BY TRANSFORMING COLORECTAL HEALTHCARE

A comprehensive overview
of services and solutions for
colorectal procedures



Medtronic
Further, Together



BETTER PATIENT OUTCOMES THROUGH LESS INVASIVE CARE

✓ INTRODUCTION

> PATIENT EXPERIENCE

> CLINICAL CHALLENGES

> PRODUCTS

We're committed to bringing a less invasive approach to Colorectal Health. One that puts the patient at the center and spans from diagnosis to treatment and recovery. And we're not going to do it alone.

We're inviting partners new and old to join us in improving the complete patient experience. Through powerful collaborations, we'll identify new, less invasive therapies — and make today's solutions even more effective.

Together, we'll bridge treatment gaps and pave the way to increased access for all patients. And deliver better care, better patient outcomes, and better lives for everyone.



Medtronic



IMPROVING THE COMPLETE PATIENT EXPERIENCE

✓ INTRODUCTION

> PATIENT EXPERIENCE

> CLINICAL CHALLENGES

> PRODUCTS



EARLY DIAGNOSIS

Early diagnosis is key to a better patient experience. Our products span the care continuum to help give your patients a less invasive journey from diagnosis to treatment and recovery.



INNOVATIVE TREATMENT

Helping you overcome challenges requires innovation.

Together, we can develop transformative technologies and create services that rise to a common goal: better patient outcomes.



FASTER RECOVERY

Shorter hospital stays, faster recoveries, and a better quality of life for the long-term. Positive patient outcomes are at the heart of everything we do.



Medtronic



LESS INVASIVE. **MORE BENEFITS.**

A shift in surgical technique is making a huge impact for patients and healthcare providers

Minimally invasive surgery (MIS) offers multiple advantages over open surgery, specifically¹⁻³:

- Fewer readmissions
- Less postoperative pain
- Lower costs

Taken together, the benefits of MIS over open surgery are significant. That's why we're working with partners to overcome infrastructure and training barriers — to make MIS the first choice in colorectal healthcare.



ADVANTAGES WITH MIS



✓ INTRODUCTION

> PATIENT EXPERIENCE

> CLINICAL CHALLENGES

> PRODUCTS



REFERENCES 

Medtronic



PATIENT FIRST SOLUTIONS FROM DIAGNOSIS TO RECOVERY

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

> OPEN STEPS



Medtronic



MIS PROCEDURAL STEPS: LESS INVASIVE TECHNOLOGY. BETTER PATIENT OUTCOMES.

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS




> MOBILIZATION & DISSECTION

> DIVISION OF CRITICAL STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS

				
ACCESS	MOBILIZATION & DISSECTION	DIVISION OF CRITICAL STRUCTURES	ANASTOMOSIS	CLOSURE
VersaStep™ Plus radially expanding trocar >	LigaSure™ retractable L-hook laparoscopic sealer/divider >	Signia™ stapling system >	EEA™ circular stapler with Tri-Staple™ technology >	Polysorb™ suture >
VersaOne™ optical trocars >	Sonicision™ cordless ultrasonic dissection device >	Endo GIA™ reloads with Tri-Staple™ technology >	EEA™ stapler with DST Series™ technology >	Biosyn™ suture >
VersaOne™ fascial closure system >	LigaSure™ Maryland jaw device >	Endo GIA™ radial reload with Tri-Staple™ technology >		Maxon™ suture >
SurgiSleeve™ wound protector >		Endo GIA™ ultra universal stapler >		V-Loc™ knotless wound closure device >



VersaStep™ Plus Radially Expanding Trocar

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

✓ ACCESS

> VersaStep™ Plus
radially expanding trocar

> VersaOne™
optical trocar

> VersaOne™ fascial
closure system

> SurgiSleeve™
wound protector

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS

Smooth insertion and secure fixation

A bladeless access device that combines the advantages of Step™ radial dilation technology with the Versaseal™ Plus universal seal.



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic





VersaOne™ Optical Trocar

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

✓ ACCESS

> VersaStep™ Plus
radially expanding trocar

> **VersaOne™
optical trocar**

> VersaOne™ fascial
closure system

> SurgiSleeve™
wound protector

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS



Access without the excess

A universal trocar system
that offers visualization on the
first entry of the tissue layer.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic





VersaOne™ Fascial Closure System

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

✓ ACCESS

> VersaStep™ Plus
radially expanding trocar

> VersaOne™
optical trocar

> **VersaOne™ fascial
closure system**

> SurgiSleeve™
wound protector

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS

A trocar. A closure system.
An all-in-one device.¹

The VersaOne™ fascial closure system is a unique all-in-one¹ solution that serves as a trocar and a fascial closure device to deliver:

- Consistent port-site closure¹
- Added procedural efficiency^{1,†}
- Ease of use¹

†Saves time compared to the competitive closure device, suture passers, and hand sutures.



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic





SurgiSleeve™ Wound Protector

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

✓ ACCESS

> VersaStep™ Plus
radially expanding trocar

> VersaOne™
optical trocar

> VersaOne™ fascial
closure system

> **SurgiSleeve™
wound protector**

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS

**Maximized exposure. Wound protection.
Superior strength.¹**

SurgiSleeve™ wound protector provides circumferential elastic retraction that maximizes working area for optimal exposure and visualization.

In a randomized clinical trial, the use of barrier wound protection in open elective colorectal resectional surgery resulted in a clinically significant reduction in incisional surgical site infections. The study concluded that barrier wound protection of this nature should be considered routine in this type of surgery.²



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



LigaSure™ Retractable L-Hook Laparoscopic Sealer/Divider

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

✓ MOBILIZATION & DISSECTION

> LigaSure™ retractable L-hook laparoscopic sealer/divider

> Sonicision™ cordless ultrasonic dissector

> LigaSure™ Maryland jaw device

> DIVISION OF CRITICAL STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS



**One instrument. Five functions.
Countless benefits.**

The functionality of five laparoscopic devices come together in one versatile instrument^{1,†} — with the reliable performance of LigaSure™ technology and Valleylab™ monopolar energy.

†29 out of 29 surgeons evaluated agreed.

TECHNOLOGY	PERFORMANCE	SPECIFICATIONS
>	>	>

REFERENCES >



Medtronic



Sonicision™

Cordless Ultrasonic Dissector

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

✓ MOBILIZATION & DISSECTION

> LigaSure™ retractable L-hook laparoscopic sealer/divider

> **Sonicision™ cordless ultrasonic dissector**


> LigaSure™ Maryland jaw device

> DIVISION OF CRITICAL STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS



Energy whenever and wherever you want it

The first cordless ultrasonic dissection system meets all your ultrasonic dissection needs and — with its cordless design — contributes to a safer operating room (OR).¹

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



LigaSure™ Maryland Jaw Device

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

✓ MOBILIZATION &
DISSECTION

> LigaSure™ retractable
L-hook laparoscopic
sealer/divider

> Sonicision™ cordless
ultrasonic dissector

> LigaSure™ Maryland
jaw device

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS



Efficient, versatile, and multifunctional

Driven by reliable LigaSure™ technology, this multifunctional tool brings together the following features in a single device¹:

- One-step sealing
- A Maryland dissector
- An atraumatic grasper
- Cold scissors

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



Signia™ Stapling System

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION &
DISSECTION

✓ DIVISION OF CRITICAL
STRUCTURES

> Signia™
stapling system

> Endo GIA™ reloads with
Tri-Staple™ technology

> Endo GIA™ radial
reload with Tri-Staple™
technology

> Endo GIA™
Ultra Universal stapler

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS



Not just a new stapler.
A smart stapler.

Smart technology that gives you
real-time feedback and powered
rotation, articulation, and firing —
with one hand.¹

TECHNOLOGY	PERFORMANCE	SPECIFICATIONS
>	>	>

REFERENCES >



Medtronic



Endo GIA™ Reloads with Tri-Staple™ Technology

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION &
DISSECTION

✓ DIVISION OF CRITICAL
STRUCTURES

> Signia™
stapling system

> Endo GIA™ reloads with
Tri-Staple™ technology

> Endo GIA™ radial
reload with Tri-Staple™
technology

> Endo GIA™
Ultra Universal stapler

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS



Stapling made smarter

Tri-Staple™ technology reloads are designed to meet critical surgical challenges.

TECHNOLOGY	PERFORMANCE	SPECIFICATIONS
>	>	>

REFERENCES >

Medtronic





Endo GIA™

Radial Reload with Tri-Staple™ Technology



> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION &
DISSECTION

✓ DIVISION OF CRITICAL
STRUCTURES

> Signia™
stapling system

> Endo GIA™ reloads with
Tri-Staple™ technology

> Endo GIA™ radial
reload with Tri-Staple™
technology

> Endo GIA™
Ultra Universal stapler

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS

Below and beyond in lower anterior resection (LAR)

Setting new standards for access,
maneuverability, and staple line security in
LAR procedures with this innovative design.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



Endo GIA™ Ultra Universal Stapler

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION &
DISSECTION

✓ DIVISION OF CRITICAL
STRUCTURES

> Signia™
stapling system

> Endo GIA™ reloads with
Tri-Staple™ technology

> Endo GIA™ radial
reload with Tri-Staple™
technology

> Endo GIA™
Ultra Universal stapler

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS



The Endo GIA™ Ultra Universal stapler combines an ergonomic design, articulation, and one-handed grasping.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



Medtronic



EEA™ Circular Stapler with Tri-Staple™ Technology



> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

✓ ANASTOMOSIS

> EEA™ circular stapler
with Tri-Staple™
technology

> EEA™ stapler with DST
Series™ technology

> CLOSURE

> OPEN STEPS

Staple line security. Times three.

The proven performance of Tri-Staple™
technology on the EEA™ circular stapler.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



EEA™ Stapler with DST Series™ Technology



> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

✓ ANASTOMOSIS

> EEA™ circular stapler
with Tri-Staple™
technology

> EEA™ stapler with DST
Series™ technology

> CLOSURE

> OPEN STEPS

The EEA™ stapler with DST Series™ technology provides reliable and consistent performance backed by years of clinical experience.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



Polysorb™

Braided Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION & DISSECTION

> DIVISION OF CRITICAL STRUCTURES

> ANASTOMOSIS

✓ CLOSURE

> Polysorb™ suture

> Biosyn™ suture

> Maxon™ suture

> V-Loc™ knotless wound closure device

> OPEN STEPS

Strength and security when you need it most

Braided synthetic absorbable suture with applications in all surgical specialities where a strong absorbable suture is needed.



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic





Biosyn™

Monofilament Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION & DISSECTION

> DIVISION OF CRITICAL STRUCTURES

> ANASTOMOSIS

✓ CLOSURE

> Polysorb™ suture

> **Biosyn™ suture**

> Maxon™ suture

> V-Loc™ knotless wound closure device

> OPEN STEPS

Our customers have trusted Biosyn™ sutures for over 20 years

Biosyn™ sutures provide the strength retention and absorption profiles to meet your surgical and tissue healing needs.



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic





Maxon™ Monofilament Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

✓ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device

> OPEN STEPS

Our customers have trusted
Maxon™ sutures for over 30 years

Maxon™ sutures are long-term
monofilament absorbable sutures
that provide excellent strength over
the critical wound healing period.



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic





V-Loc™

Knotless Wound Closure Device

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION & DISSECTION

> DIVISION OF CRITICAL STRUCTURES

> ANASTOMOSIS

✓ CLOSURE

> Polysorb™ suture

> Biosyn™ suture


> Maxon™ suture

> V-Loc™ knotless wound closure device

> OPEN STEPS

The V-Loc™ knotless wound closure device is a revolutionary technology that closes wounds securely without the need to tie knots

The device achieves a secure closure for patients.



TECHNOLOGY	PERFORMANCE	SPECIFICATIONS
>	>	>



Medtronic



OPEN PROCEDURAL STEPS: TECHNOLOGY YOU CAN TRUST.

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS





✓ OPEN STEPS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

			
MOBILIZATION & DISSECTION	DIVISION OF CRITICAL STRUCTURES	ANASTOMOSIS	CLOSURE
LigaSure™ Maryland jaw device >	TA™ stapler with DST Series™ technology >	EEA™ circular stapler with Tri-Staple™ technology >	Polysorb™ suture >
LigaSure Impact™ curved large jaw open sealer-divider, nano-coated >	GIA™ stapler with DST Series™ technology >	EEA™ stapler with DST Series™ technology >	Biosyn™ suture >
Sonicision™ cordless ultrasonic dissection device >	SurgiSleeve™ wound protector >		Maxon™ suture >
			V-Loc™ knotless wound closure device >



LigaSure™ Maryland Jaw Device

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

✓ MOBILIZATION &
DISSECTION

> LigaSure™ Maryland
jaw device

> LigaSure Impact™
open instrument

> Sonicision™ cordless
ultrasonic dissector

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE



Efficient, versatile, and multifunctional

Driven by reliable LigaSure™ technology, this multifunctional tool brings together the following features in a single device¹:

- One-step sealing
- A Maryland dissector
- An atraumatic grasper
- Cold scissors

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



LigaSure Impact™

Curved Large Jaw Sealer/ Divider, Nano-coated

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

✓ MOBILIZATION &
DISSECTION

> LigaSure™ Maryland
jaw device

> LigaSure Impact™
open instrument

> Sonicision™ cordless
ultrasonic dissector

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE



Seal smarter. Seal better.

Multifunctional tissue fusion system with integrated cutting mechanism — delivers reliable performance and procedural efficiency.¹

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



Sonicision™

Cordless Ultrasonic Dissector

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

✓ MOBILIZATION & DISSECTION

> LigaSure™ Maryland jaw device

> LigaSure Impact™ open instrument

> **Sonicision™ cordless ultrasonic dissector**

> DIVISION OF CRITICAL STRUCTURES

> ANASTOMOSIS

> CLOSURE



Energy whenever and wherever you want it

The first cordless ultrasonic dissection system meets all your ultrasonic dissection needs and — with its cordless design — contributes to a safer OR.¹

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



TA™ Stapler with DST Series™ Technology

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION &
DISSECTION

✓ DIVISION OF CRITICAL
STRUCTURES


> TA™ stapler with DST
Series™ technology

> GIA™ stapler with DST
Series™ technology

> SurgiSleeve™
wound protector

> ANASTOMOSIS

> CLOSURE



Reliability in staple formation with
our proprietary directional staple
technology (DST)

A reloadable linear stapler that places
one double staggered row of titanium
staples to occlude tissue during open
colorectal surgery.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



GIA™ Stapler with DST Series™ Technology

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION &
DISSECTION

✓ DIVISION OF CRITICAL
STRUCTURES

> TA™ stapler with DST
Series™ technology

> GIA™ stapler with DST
Series™ technology

> SurgiSleeve™
wound protector

> ANASTOMOSIS

> CLOSURE



Reliability in staple formation with our proprietary DST

A reloadable linear stapler that places two double-staggered rows of titanium staples and simultaneously cuts and divides tissue during open colorectal surgery.

TECHNOLOGY	PERFORMANCE	SPECIFICATIONS
>	>	>

REFERENCES >



Medtronic



SurgiSleeve™ Wound Protector

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION &
DISSECTION

✓ DIVISION OF CRITICAL
STRUCTURES

> TA™ stapler with DST
Series™ technology

> GIA™ stapler with DST
Series™ technology

> SurgiSleeve™
wound protector

> ANASTOMOSIS

> CLOSURE

**Maximized exposure. Wound protection.
Superior strength.¹**

SurgiSleeve™ wound protector provides circumferential elastic retraction that maximizes working area for optimal exposure and visualization.

In a randomized clinical trial, the use of barrier wound protection in open elective colorectal resectional surgery resulted in a clinically significant reduction in incisional surgical site infections. The study concluded that barrier wound protection of this nature should be considered routine in this type of surgery.²



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



EEA™ Circular Stapler with Tri-Staple™ Technology



> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

✓ ANASTOMOSIS

> EEA™ circular stapler
with Tri-Staple™
technology

> EEA™ stapler with DST
Series™ technology

> CLOSURE

Staple line security. Times three.

The proven performance of Tri-Staple™
technology on the EEA™ circular stapler.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



EEA™ Stapler with DST Series™ Technology



> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

✓ ANASTOMOSIS

> EEA™ circular stapler
with Tri-Staple™
technology

> EEA™ stapler with DST
Series™ technology

> CLOSURE

The EEA™ stapler with DST Series™ technology provides reliable and consistent performance backed by years of clinical experience.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



Polysorb™ Braided Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

✓ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device

Strength and security when you need it most

Braided synthetic absorbable
suture with applications in all
surgical specialities where a strong
absorbable suture is needed.



TECHNOLOGY	PERFORMANCE	SPECIFICATIONS
>	>	>

REFERENCES >

Medtronic





Biosyn™ Monofilament Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

✓ CLOSURE

> Polysorb™
suture

> **Biosyn™**
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device

Our customers have trusted
Biosyn™ sutures for over 20 years

Biosyn™ sutures provide the
strength retention and absorption
profiles to meet your surgical and
tissue healing needs.



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



Maxon™ Monofilament Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

✓ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device

Our customers have trusted
Maxon™ sutures for over 30 years

Maxon™ sutures are long-term
monofilament absorbable sutures
that provide excellent strength over
the critical wound healing period.



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



V-Loc™

Knotless Wound Closure Device

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION & DISSECTION

> DIVISION OF CRITICAL STRUCTURES

> ANASTOMOSIS

✓ CLOSURE

> Polysorb™ suture


> Biosyn™ suture

> Maxon™ suture

> V-Loc™ knotless wound closure device

The V-Loc™ knotless wound closure device is a revolutionary technology that closes wounds securely without the need to tie knots

The device achieves a secure closure for patients.



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



Medtronic



LESS INVASIVE. MORE BENEFITS.

A shift in surgical technique is making a huge impact for patients and healthcare providers

INTRODUCTION

PATIENT EXPERIENCE



SHORTER HOSPITAL STAYS.
2.5 DAYS DIFFERENCE, ON AVERAGE.^{4-8, 11-19}



53% LESS BLOOD LOSS,
ON AVERAGE.⁵⁻¹⁰



FASTER RECOVERY.
QUICKER RETURN TO NORMAL ACTIVITY.²⁰



3.2% LOWER RATE
OF SURGICAL SITE INFECTION.²¹⁻²³



REFERENCES

Medtronic



REFERENCES

1. Fitch K, Engel T, Bochner A. Cost differences between open and minimally invasive surgery. *Managed Care*. 2015;24(9):40–48.
2. Tiwari MM, Reynoso JF, High R, Tsang AW, Oleynikov D. Safety, efficacy, and cost effectiveness of common laparoscopic procedures. *Surg Endosc*. 2011;25(4):1127–1135.
3. Roumm AR, Pizzi L, Goldfarb NI, Cohn H. Minimally invasive: minimally reimbursed? An examination of six laparoscopic surgical procedures. *Surg Innovation*. 2005;12(3):261–287.
4. Juo YY, Hyder O, Haider AH, Camp M, Lidor A, Ahuja N. Is minimally invasive colon resection better than traditional approaches?: First comprehensive national examination with propensity score matching. *JAMA Surg*. 2014;149(2):177–184.
5. Veldkamp R, Kuhry E, Hop WC, et al. Laparoscopic surgery versus open surgery for colon cancer: short-term outcomes of a randomised trial. *Lancet Oncol*. 2005;6(7):477–484.
6. Orcutt ST, Marshall CL, Robinson CN, et al. Minimally invasive surgery in colon cancer patients leads to improved short-term outcomes and excellent oncologic results. *Am J Surg*. 2011;202(5):528–531.
7. Chen K, Zhang Z, Zuo Y, Ren S. Comparison of the clinical outcomes of laparoscopic-assisted versus open surgery for colorectal cancer. *Oncol Lett*. 2014;7(4):1213–1218.
8. da Luz Moreira A, Kiran RP, Kirat HT, et al. Laparoscopic versus open colectomy for patients with American Society of Anesthesiology (ASA) classifications 3 and 4: the minimally invasive approach is associated with significantly quicker recovery and reduced costs. *Surg Endosc*. 2010;24(6):1280–1286.
9. Braga M, Frasson M, Zuliani W, Vignali A, Pecorelli N, Di Carlo V. Randomized clinical trial of laparoscopic versus open left colonic resection. *Br J Surg*. 2010;97(8):1180–1186.
10. Li Z, Li D, Jie Z, Zhang G, Liu Y. Comparative study on therapeutic efficacy between hand-assisted laparoscopic surgery and conventional laparotomy for acute obstructive right-sided colon cancer. *J Laparoendosc Adv Surg Tech A*. 2015;25(7):548–554.
11. Steele SR, Brown TA, Rush RM, Martin MJ. Laparoscopic vs. open colectomy for colon cancer: results from a large nationwide population-based analysis. *J Gastrointest Surg*. 2008;12(3):583–591.





REFERENCES

12. Odermatt M, Miskovic D, Siddiqi N, Khan J, Parvaiz A. Short- and long-term outcomes after laparoscopic versus open emergency resection for colon cancer: an observational propensity score-matched study. *World J Surg.* 2013;37(10):2458–2467.
13. Prakash K, Varma D, Rajan M, et al. Laparoscopic colonic resection for rectosigmoid colonic tumours: a retrospective analysis and comparison with open resection. *Indian J Surg.* 2010;72(4):318–322.
14. Hardy KM, Kwong J, Pitzul KB, et al. A cost comparison of laparoscopic and open colon surgery in a publicly funded academic institution. *Surg Endosc.* 2014;28(4):1213–1222.
15. Thompson BS, Coory MD, Gordon LG, Lumley JW. Cost savings for elective laparoscopic resection compared with open resection for colorectal cancer in a region of high uptake. *Surg Endosc.* 2014;28(5):1515–1521.
16. Marshall CL, Chen GJ, Robinson CN, et al. Establishment of a minimally invasive surgery program leads to decreased inpatient cost of care in veterans with colon cancer. *Am J Surg.* 2010;200(5):632–635.
17. Day AR, Smith RV, Jourdan IC, Rockall TA. Survival following laparoscopic and open colorectal surgery. *Surg Endosc.* 2013;27(7):2415–2421.
18. Kapritsou M, Korkolis DP, Konstantinou EA. Open or laparoscopic surgery for colorectal cancer: a retrospective comparative study. *Gastroenterol Nurs.* 2013;36(1):37–41.
19. Agarwal S, Gincherman M, Birnbaum E, Fleshman JW, Mutch M. Comparison of long-term follow up of laparoscopic versus open colectomy for transverse colon cancer. *Proc (Bayl Univ Med Cent).* 2015;28(3):296–299.
20. Society of American Gastrointestinal and Endoscopic Surgeons (SAGES) Webmaster. Laparoscopic colon resection surgery patient information from SAGES website. <https://www.fascrs.org/patients/disease-condition/minimally-invasive-surgery-expanded-version>. Updated March 1, 2015. Accessed January 10, 2017.
21. Bilimoria KY, Bentrem DJ, Merkow RP, et al. Laparoscopic-assisted vs. open colectomy for cancer: comparison of short-term outcomes from 121 hospitals. *J Gastrointest Surg.* 2008;12(11):2001–2009.
22. Kiran RP, Kirat HT, Ozturk E, Geisler DP, Remzi FH. Does the learning curve during laparoscopic colectomy adversely affect costs? *Surg Endosc.* 2010;24(11):2718–2722.
23. Wilson MZ, Hollenbeak CS, Stewart DB. Laparoscopic colectomy is associated with a lower incidence of postoperative complications than open colectomy: a propensity score-matched cohort analysis. *Colorectal Dis.* 2014;16(5):382–389.



VersaStep™ Plus Radially Expanding Trocar

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

✓ ACCESSORIES

> VersaStep™
radially expanding

> VersaStep™
optical trocar

> VersaStep™
closure trocar

> SurgiSeal™
wound protector

> MOBILIZATION
DISSECTION

> DIVISION
STRUCTURE

> ANASTOMOSES

> CLOSURE

> OPEN STEPS

TECHNOLOGY

Reliability every step of the way:

- Step™ radially expanding technology allows the VersaStep™ trocar to yield smaller fascial defects for an equivalent cannula size — compared to conventional bladed trocars.¹
- A durable seal ensures that insufflation is maintained even during long advanced cases needing many instrument exchanges.
- A rotating instrument seal provides excellent tactile feel and full range of motion while maintaining an airtight seal.



**NO
CONVERTER
ACCESSORIES
ARE REQUIRED
FOR
INSTRUMENTS
OF DIFFERENT
SIZES**



REFERENCES



Medtronic



PERFORMANCE

CHALLENGE

SOLUTION

Risk of tissue trauma.

Step™ radially expanding technology allows the VersaStep™ trocar to yield smaller fascial defects for an equivalent cannula size — compared to conventional bladed trocars.¹

Need for converter accessories for instruments of different sizes.

VersaStep™ Plus radially expanding trocar maintains pneumoperitoneum while facilitating exchange of instruments ranging from 4.5 mm to 12 mm.

Lack of flexibility during instrument exchange.

VersaStep™ Plus self adjusts to accept 5 mm to 12 mm instruments minimizing the loss of safe working space.

VersaStep™ Plus Radially Expanding Trocar

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

∨ ACCESS

> VersaStep™ Plus
radially expanding trocar

> VersaOne™
optical trocar

> VersaOne™ fascial
closure system

> SurgiSleeve™
wound protector

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS



SPECIFICATIONS



Available in 5 mm, 11 mm,
12 mm, and 15 mm sizes.

A diaphragm access device that
combines the advantages of
Step™ radial dilation technology
with the Versaseal™ Plus
universal seal.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



VersaStep™ Plus Radially Expanding Trocar

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

∨ ACCESS

> VersaStep™ Plus
radially expanding tro

> VersaOne™
optical trocar

> VersaOne™ fascial
closure system

> SurgiSleeve™
wound protector

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS

REFERENCE

1. Johnson WH, Fecher AM, McMahon RL, McMahon RL, Grant JP, Pryor AD. VersaStep™ trocar hernia rate in unclosed fascial defects in bariatric patients. *Surg Endosc.* 2006;20(10):1584–1586.

combines the advantages of
Step™ radial dilation technology
with the Versaseal™ Plus
universal seal.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



VersaOne™ Optical Trocar

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

∨ ACCESS

> VersaStep
radially exp

> VersaOne
optical tro

> VersaOne
closure sys

> SurgiSleeve
wound pro

> MOBILIZATI
DISSECTIO

> DIVISION O
STRUCTURE

> ANASTOMO

> CLOSURE

> OPEN STEPS

TECHNOLOGY

VersaOne™ optical trocar delivers:

- "Dolphin-nose tip" to facilitate smoother insertion for easier trocar placement versus Applied Medical¹
- Advanced fixation ribs to keep cannula from moving during surgical procedure¹
- Excellent visualization during insertion



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

PERFORMANCE



CHALLENGE

SOLUTION

Risk of injury during primary port placement.

VersaOne™ optical trocar has increased visibility with optical tip geometry allowing for an informed entry.^{1,†}

Insufficient visualization.

VersaOne™ optical trocar has a clear cannula for visualization.[†]

Difficult insertion.

VersaOne™ optical trocar has a bladeless "dolphin nose" obturator tip that facilitates smooth insertion for easier trocar placement, as well as a bladed obturator option.¹

Multiple codes and SKUs across procedures.

VersaOne™ optical trocar has an interchangeable obturator allowing flexibility for procedural needs.

†Images recorded on March 14, 2013, as part of a trocar clarity evaluation in an inanimate abdominal model; images representative of evaluation results.

VersaOne™ Optical Trocar

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

∨ ACCESS

> VersaStep™ Plus
radially expanding trocar

> VersaOne™
optical trocar

> VersaOne™ fascial
closure system

> SurgiSleeve™
wound protector


> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS



SPECIFICATIONS



Available in 5 mm, 11 mm,
12 mm, and 15 mm sizes.

Access without the excess

A universal trocar system
that offers visualization on the
first entry of the tissue layer.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

VersaOne™ Optical Trocar

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

∨ ACCESS

> VersaStep™ Plus
radially expanding tro

> VersaOne™
optical trocar

> VersaOne™ fascial
closure system

> SurgiSleeve™
wound protector

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS

REFERENCES

1. Based on internal test report #2143-114, Applies to Medtronic 12 mm trocars, when compared to Applied Kii™* (12 mm trocar, Z-thread cannula). March 2013.

A universal trocar system
that offers visualization on the
first entry of the tissue layer.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

VersaOne™ Fascial Closure System

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

∨ ACCESS

> VersaStep
radially exp

> VersaOne™
optical troc

> VersaOne™
closure sys

> SurgiSleeve
wound pro

> MOBILIZATI
DISSECTIO

> DIVISION O
STRUCTURE

> ANASTOMO

> CLOSURE

> OPEN STEPS

TECHNOLOGY

Access and closure with a single device

The VersaOne™ fascial closure system is an all-in-one¹ device that brings consistency to laparoscopic closures.¹ It also:

- Eliminates the need to remove the trocar before closing^{1,2}
- Makes reinsufflation unnecessary — pneumoperitoneum can be maintained throughout the procedure²



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

PERFORMANCE



CHALLENGE

SOLUTION

Risk of injury during primary port placement.

VersaOne™ fascial closure system has increased visibility with optical tip geometry allowing for an informed entry.^{3,†}

Insufficient visualization.

VersaOne™ fascial closure system has a clear cannula for visualization.[†]

Difficult insertion.

VersaOne™ fascial closure system has a bladed and bladeless "dolphin-nose" obturator tip facilitating smooth insertion for easier trocar placement.^{4,†}

Multiple codes and SKUs needs across procedures.

VersaOne™ fascial closure system has an interchangeable obturator allowing flexibility for procedural needs.

†As compared to our legacy trocars.

VersaOne™ Fascial Closure System

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

∨ ACCESS

> VersaStep™ Plus
radially expanding trocar

> VersaOne™
optical trocar

> VersaOne™ fascial
closure system

> SurgiSleeve™
wound protector

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS

A trocar. A closure system.
An all-in-one device.¹

SPECIFICATIONS

Available in 12 mm optical and
bladeless models.

- Added procedural efficiency^{1,†}
- Ease of use¹

†Saves time compared to the competitive closure
device, suture passers, and hand sutures.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

VersaOne™ Fascial Closure System

> INTRODUCTION

∨ PRODUCT

∨ MISSTEP

∨ ACCESS

>

>

>

>

> MOE
DISS

> DIVI
STR

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS

REFERENCES

1. Based on internal validation report #RE00098009 Rev A, DoubleTime VersaOne™ fascial closure system Miami lab VOC report: 13 out of 16 surgeons surveyed agreed. March 28, 2017.
2. Based on internal validation report #RE00062341 Rev A, Summative usability of project DoubleTime VersaOne™ fascial closure system. August 24–25, 2016.
3. Images recorded on March 14, 2013 as part of a trocar clarity evaluation in an inanimate abdominal model; images representative of evaluation results.
4. Based on internal test report #2143-114, Applies to Medtronic 12 mm trocars, when compared to Applied Kii™* (12 mm trocar, Z-thread cannula). March 2013.

†Saves time compared to the competitive closure device, suture passers, and hand sutures.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

TECHNOLOGY

SurgiSleeve™ wound protector provides optimal exposure and wound protection in surgical cases

SurgiSleeve™ wound protectors:

- Are superior in film strength to Alexis™* wound protector/retractors¹
- Offer equivalent or greater exposure versus Alexis™* O with a retraction ring in the large and extra-large sizes as rated by 98 percent of surgeons surveyed^{3,4}
- Are available with large and extra-large retraction rings that are fast to place and remove⁵
- Are faster to roll down compared to Alexis™* O as rated by a majority of surgeons surveyed⁶
- Come in 5 SKUs to cover all your wound protection needs versus 15 SKUs for Alexis™* O wound protector/retractor





PERFORMANCE

CHALLENGE

SOLUTION

Sufficient retraction at incision site.^{5,7-10}

SurgiSleeve™ wound protector has a blue proximal ring that is not only easier to roll down compared to Alexis™* O, but can be rolled down by just one person.†

Adequate exposure to the body.^{5,7-10}

SurgiSleeve™ wound protector has a retraction ring (large and extra-large sizes) providing optimal retraction and visibility compared to Alexis™* O and a stable platform for access to the body.⁶

Potential for wound contamination.^{5,7-10}

SurgiSleeve™ wound protector has three times stronger film material¹ compared to Alexis™* O offering optimal wound protection.

†Survey of 29 General and OB/GYN surgeons of SurgiSleeve™ Wound Protector and Alexis®* O conducted during a product demonstration event October 28–29, 2014 in San Francisco, CA.

SurgiSleeve™ Wound Protector

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

∨ ACCESS

> VersaStep™ Plus
radially expanding trocar

> VersaOne™
optical trocar

> VersaOne™ fascial
closure system

> SurgiSleeve™
wound protector

> MOBILIZATION &
DISSECTION


> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS

Maximized exposure. Wound protection.
Superior strength.¹



SPECIFICATIONS

Available in 5 sizes: extra-small, small, medium, large with retraction ring, and extra-large with retraction ring.

clinically significant reduction in incisional surgical site infections. The study concluded that barrier wound protection of this nature should be considered routine in this type of surgery.²

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



REFERENCES

1. Based on internal report #2151-002, Comparison of mean film material strengths between SurgiSleeve™ wound protector and Alexis™* wound protector, as part of a puncture resistance material evaluation conducted on February 9, 2012. Statistically significant results $p = 0.004$.
2. Reid K, Pockney P, Draganic B, Smith SR. Barrier wound protection decreases surgical site infection in open elective colorectal surgery: A randomized clinical trial. *Dis Colon Rectum*. 2010; 53: 1374–1380.
3. Based on internal report #2151-021, SurgiSleeve™ large wound protector with retraction ring equivalence testing. October 30, 2014, and
4. Based on internal report #2151-023, SurgiSleeve™ extra-large wound protector with retraction ring equivalence testing. November 4, 2014.
5. Voice-of-the-customer survey of 29 general and OB/GYN surgeons of SurgiSleeve™ wound protector and Alexis™* O large and extra-large conducted during a product demonstration event, October 28–29, 2014, in San Francisco, CA.
6. Information in the table was based on (1) Survey of 29 General and OB/GYN surgeons of SurgiSleeve™ wound protector and Alexis™* O large and extra-large conducted during a product demonstration event Oct. 28–29, 2014 in San Francisco, CA; and (2) the animal labs for the large & extra-large ring with five general surgeons and three OB/GYN surgeons on Aug. 5–6 in the North Haven, CT, USA facility.
7. Edwards JP, Ho AL, Tee MC, Dixon E, Ball CG. Wound protectors reduce surgical site infection: a meta-analysis of randomized controlled trials. *Ann Surg*. 2012; 256(1): 53–59.
8. Horiuchi T, Tanishima H, Tamagawa K, et al. A wound protector shields incision sites from bacterial invasion. *Surg Infect (Larchmt)*. 2010;11(6):501–503.
9. Horiuchi T, Tanishima H, Tamagawa K, et al. Randomized, controlled investigation of the anti-infective properties of the Alexis retractor/protector of incision sites. *J Trauma*. 2007;62(1):212–215.
10. Lee P, Waxman K, Taylor B, Yim S. Use of wound-protection system and postoperative wound-infection rates in open appendectomy. *Arch Surg*. 2009;144(9): 872–875.


LigaSure™ Retractable L-Hook Laparoscopic Sealer/Divider

TECHNOLOGY

The only device of its kind, it delivers the benefits of²:

- Precise Valleylab™ monopolar dissection^{1,†}
- One-step LigaSure™ vessel sealing^{1,†}
- Atraumatic grasping²
- Cold cutting²
- Maryland-style blunt dissection²

†29 out of 29 surgeons evaluated agreed.



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



PERFORMANCE

CHALLENGE

Laparoscopic colorectal procedures take time and require frequent position resets after instrument exchanges.

Colorectal cases often require access to tight spaces deep in the pelvis.

SOLUTION

The LigaSure™ retractable L-hook device provides enhanced surgical flow due to one-step activation and decreased instrument exchanges.^{1,†}

The LigaSure™ retractable L-hook device's tapered, curved jaw profile allows for access in tight spaces.^{1,†}

†29 out of 29 surgeons evaluated agreed.

LigaSure™

Retractable L-Hook Laparoscopic Sealer/Divider

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

> ACCESS

∨ MOBILIZATION &
DISSECTION

> LigaSure™ retractable
L-hook laparoscopic
sealer/divider

> Sonicision™ cordless
ultrasonic dissector

> LigaSure™ Maryland
jaw device

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS



SPECIFICATIONS



LigaSure™ technology can seal:

- Vessels up to and including 7 mm
- Lymphatics
- Tissue bundles

Available in 37 cm and 44 cm sizes.

and Valleylab™ monopolar energy.

†29 out of 29 surgeons evaluated agreed.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

LigaSure™

Retractable L-Hook Laparoscopic Sealer/Divider

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

> ACCESS

∨ MOBILIZATION &
DISSECTION

> LigaSure™ retractable
L-hook laparoscopic
sealer/divider

> Sonicision™ cordless
ultrasonic dissector

> LigaSure™ Maryland
jaw device

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS

REFERENCES

1. Based on internal test report #RE00041188, Validation marketing claims: independent surgeon feedback collected during Medtronic-sponsored cadaver and porcine labs. October 2015 and February 2016.
2. Based on internal test report #RE00032739 Rev A, LigaSure™ Hook (LF5637, LF5644) validation surgeon evaluation report: independent surgeon feedback collected during Medtronic-sponsored cadaver and porcine lab conducted in Houston, Texas. February 2016.

and Valleylab™ monopolar energy.

†29 out of 29 surgeons evaluated agreed.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Sonicision™ Cordless Ultrasonic Dissector

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

✓ MOBILIZATION
DISSECTION

> LigaSure™
L-hook laparoscopic
sealer/divider

> Sonicision™
ultrasonic dissector

> LigaSure™
jaw device

> DIVISION OF
STRUCTURE

> ANASTOMOSES

> CLOSURE

> OPEN STEPS

TECHNOLOGY

Offering improved freedom of movement and¹⁻³:

- Innovative cordless technology
- Rapid dissection of tissue and security of hemostasis for vessels up to 5 mm in diameter⁴
- Intuitive energy activation control improves focus on the procedure¹⁻³



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Sonicision™

Cordless Ultrasonic Dissector



PERFORMANCE

CHALLENGE

The growing popularity of laparoscopic and robotic procedures require additional equipment in the OR.

Capital equipment requires planning and the appropriate budgets.

SOLUTION

Sonicision™ cordless ultrasonic device eliminates a cord and it has been shown that the reduction of cords improves safety in the OR.¹

Sonicision™ cordless ultrasonic device platform is easy to adopt without a large capital outlay.

Sonicision™

Cordless Ultrasonic Dissector

> INTRODUCTION

∨ PRO

∨ M



SPECIFICATIONS

Measure	Sonicision™ system	Ethicon Harmonic ACE™†	Ethicon Harmonic ACE+™‡	Olympus Thunderbeat™§
Shaft length	13 cm	14 cm	23 cm	10 cm
	26 cm	23 cm	36 cm	20 cm
	39 cm	36 cm		35 cm
	48 cm	45 cm		45 cm
Method of activation	Hand	Hand	Hand	Hand
	1 button, 2 stages	2 buttons	2 buttons	2 buttons
Frequency	55.5 kHz	55.5 kHz	55.5 kHz	47.7 kHz
Jaw length	14.5 mm	11.1 mm	10.7 mm	16.5 mm
Total weight of device assembled (with cords and transducer)	394.25 g	337.47 g	341.86 g	386.86 g

†Ethicon Harmonic ACE™* instructions for use.
 ‡Ethicon Harmonic ACE+™* instructions for use.
 §Olympus Thunderbeat™* instructions for use.

TECHNOLOGY PERFORMANCE SPECIFICATIONS



REFERENCES >

Sonicision™

Cordless Ultrasonic Dissector

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

> ACCESS

∨ MOBILIZATION & DISSECTION

> LigaSure™ retractable L-hook laparoscopic sealer/divider

> Sonicision™ cordless ultrasonic dissector

> LigaSure™ Maryland jaw device

> DIVISION OF CRITICAL STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS

REFERENCES

1. Based on internal report #R0015050, CCR review for project Milwaukee. Nov. 21, 2013.
2. Based on internal test report #R0014853 Rev A. Usability validation report surgical evaluation of Milwaukee: independent surgeon feedback collected during internally-sponsored cadaver and porcine labs. April 2010 and May 2010.
3. Based on internal report #2-108-10, Validation of Milwaukee USG1, 5 mm cordless laparoscopic ultrasonic shears. April–May, 2010.
4. Tsirline VB, Lau KN, Swan RZ, et al. Evaluation of an innovative, cordless ultrasonic dissector. *Surg Innov.* 2013 Oct;20(5):524–529.
5. Brogmus G, Leone W, Butler L, Hernandez E. Best practices in OR suite layout and equipment choices to reduce slips, trips, and falls. *AORN J.* 2007;86(3):384–394.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS

REFERENCES

Medtronic

LigaSure™ Maryland Jaw Device

TECHNOLOGY

Reliable vessel-sealing technology:

- Nonstick nano-coated jaws
- Enhanced blunt dissection compared to straight jaws^{2,†}
- Improved tip visualization compared to straight jaws^{3,‡}

The LigaSure™ Maryland jaw device (LF1944) reduces sticking by more than 38 times compared to the Ethicon ENSEAL™* G2 device.^{4,§}

†23 of 32 surgeons surveyed agreed.

‡31 of 33 surgeons surveyed agreed.

§Tissue sticking to device jaws instances measure over 110 seals per device (ForceTriad™ energy platform).



> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

✓ MOBILIZATION
DISSECTION

> LigaSure
L-hood
sealer

> Sonic
ultrasonic

> LigaSure
jaw device

> DIVISION
STRUCTURE

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS



REFERENCES



Medtronic

LigaSure™ Maryland Jaw Device



PERFORMANCE

CHALLENGE

Appropriate skeletonization of IMA is needed to allow confidence when sealing.

Mobilization of the colon requires dissection and vessel sealing which could create the need for frequent instrument exchanges.

SOLUTION

The curved jaws of the LigaSure™ Maryland jaw device allow for easy skeletonization of vessels.^{3,†}

The LigaSure™ Maryland jaw device is a one-step vessel sealer that provides the reliability, control, and consistency² of LigaSure™ technology.

†30 out of 33 surgeons surveyed agreed.



REFERENCES

Medtronic

LigaSure™ Maryland Jaw Device

> INTRODUCTION

∨ PRODUCTS



SPECIFICATIONS



	LigaSure™ Maryland jaw device	Ethicon Harmonic ACE+™†
Seal plate length	20.3 mm	13.5 mm
Cut length	18.5 mm	13.5 mm
Jaw aperture	12.7 mm	11.5 mm
Vessel size	7 mm	5 mm
Curved jaw	Yes	Yes
Degrees of curve	22 degrees	17.4 degrees
Shaft lengths	23 cm, 37 cm, 44 cm	14 cm, 23 cm, 36 cm, 45 cm
Ability to cold cut	Yes	No

†Ethicon Harmonic ACE+™* instructions for use.

▪ Cold scissors

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

LigaSure™

Maryland Jaw Device

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

✓ MOBILIZATION & DISSECTION

> LigaSure™ retractable L-hook laparoscopic sealer/divider

> Sonicision™ cordless ultrasonic dissector

> LigaSure™ Maryland jaw device

> DIVISION OF CRITICAL STRUCTURES

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS



REFERENCES

1. Based on LigaSure™ Maryland Jaw Sealer/Divider, Nano-coated [instructions for use]. Boulder, CO: Medtronic; 2016.
2. Based on internal test report #R0035742, Maryland validation labs: Independent surgeon feedback collected during porcine labs in Houston and Los Angeles. April 16 to 18 and April 30 to May 3, 2013.
3. Based on internal test report #RE00071598. Maryland Validations Labs Marketing report. Houston and Los Angeles: independent surgeon feedback collected during porcine labs. April 16–18 and April 30–May 3, 2013.
4. Based on internal test report #RE00073194, Tissue sticking comparison of the Ethicon G2™*. Voyant™* 5 mm Fusion, LigaSure™ LF1737, and LigaSure™ LF1937 devices conducted on porcine tissue using the Force Triad™ energy platform. January 18, 2017.
5. Based on internal test report #RE00140529 rev. A, LigaSure™ Maryland device, nano-coated (LF19XX) tissue testing (memo). March 5, 2018.



REFERENCES >

Medtronic

Signia™ Stapling System

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION
DISSECTION

✓ DIVISION OF
STRUCTURE

> Signia™
stapling system

> Endo GIA™
Tri-Staple™

> Endo GIA™
reload with
technology

> Endo GIA™
Ultra Univer-

> ANASTOMO-

> CLOSURE

> OPEN STEPS

TECHNOLOGY

The Signia™ stapler:

- Delivers fully powered articulation, rotation, clamping, and firing to provide precision and maneuverability¹⁻³
- Sets firing speed based on tissue clamped^{1,4-6}
- Provides an OLED screen to display real-time feedback
- Allows for single-handed operation that frees your other hand when stapling and stay focused on the surgical site¹



REFERENCES

Medtronic



PERFORMANCE

CHALLENGE

Laparoscopic colorectal surgery is challenging — especially when trying to optimally position the stapler low in the pelvis.

Poor staple formation is a risk in dense and variable tissues, especially in the colon and rectum following radiation therapy.

It's difficult to know whether you have selected the right staple reload when visualization of the tissue is limited deep in the pelvis.

It can be difficult to adequately manage multiple instruments in a laparoscopic procedure.

SOLUTION

Fully powered articulation, rotation, clamping, and firing provides greater precision and maneuverability compared to manual staplers.¹⁻³

The Signia™ stapler adjusts firing speed based on force feedback during clamp and firing to optimize stapling formation and deliver consistent staple lines.^{1,4-6}

The Signia™ stapler displays real-time feedback on an easy-to-understand OLED screen.^{1,7}

The Signia™ stapler frees your other hand to focus on the surgical site.¹

Signia™ Stapling System

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

> ACCESS

> MOBILIZATION &
DISSECTION

∨ DIVISION OF CRITICAL
STRUCTURES

> Signia™
stapling system

> Endo GIA™ reloads with
Tri-Staple™ technology

> Endo GIA™ radial
reload with Tri-Staple™
technology

> Endo GIA™
Ultra Universal stapler

> ANASTOMOSIS

> CLOSURE

> OPEN STEPS



SPECIFICATIONS



View full specifications [here](#)

Smart technology that gives you real-time feedback and powered rotation, articulation, and firing — with one hand.¹

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



Signia™ Stapling System

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

> ACCESS

> MOBILIZATION
DISSECTION

∨ DIVISION OF
STRUCTURE

> Signia™
stapling system

> Endo GIA™
Tri-Staple

> Endo GIA™
reload with
technology

> Endo GIA™
Ultra Univer

> ANASTOMOSES

> CLOSURE

> OPEN STEPS



REFERENCES

1. Based on internal test report #RE00024826, Signia™ stapling system summative usability report, Rev A. January 2016.
2. Based on internal test report #R0033554, iDrive™ Ultra PCR summary report. November 2012.
3. Based on internal test report #RE00022065, UCONN biodynamics final report on results focusing on biomechanical exposure related to laparoscopic stapler use. 2012.
4. Based on internal engineering report #PCG022, Comparison of the iDrive™-S, iDrive™-F, Ethicon 60 and universal 60 staple formation in thick tissue stomach firings. January 2016.
5. Based on internal test report #R2146-173-0, ASA verification testing with slow speed force limit evaluation. 2015.
6. Based on internal test report #R2146-151-0, Powered stapling firing speed DOE analysis and ASA parameters. 2015.
7. PT00002451 Signia™ Stapler User Manual. Page 13. Minneapolis, MN: Medtronic, 2016.



REFERENCES



Medtronic

Endo GIA™ Reloads with Tri-Staple™

TECHNOLOGY

Preserving tissue integrity is the ultimate premise behind the advanced Tri-Staple™ technology design. With its stepped cartridge face, fixed anvil, and varied height staples, Tri-Staple™ technology:

- Puts less stress on tissue during compression and clamping^{1,†}
- Allows greater perfusion into the staple line^{2,†}
- Provides superior performance in variable tissue thicknesses^{3,†}

†Preclinical results may not correlate with clinical performance in humans.



LESS STRESS

on tissue during compression and clamping^{1,†}



GREATER PERFUSION

into the staple line^{2,†}



SUPERIOR PERFORMANCE

in variable tissue thicknesses^{3,†}

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION
DISSECTION

✓ DIVISION OF
STRUCTURE

> Signia™
stapling system

> Endo GIA™
Tri-Staple™

> Endo GIA™
reload with
technology

> Endo GIA™
Ultra Univer-

> ANASTOMO-

> CLOSURE

> OPEN STEPS



REFERENCES



Medtronic



PERFORMANCE

CHALLENGE

SOLUTION

Ischemia and reduced perfusion may contribute to anastomotic leaks.⁷

The profile of Tri-Staple™ technology design allows greater perfusion into the staple line when compared to conventional stapling technologies.^{2,†}

Maintaining hemostasis while minimizing risk of ischemia during colorectal surgery is critical to reducing intraoperative blood loss.

Endo GIA™ reloads have a unique stepped cartridge face and their varied height staples approximate tissue to the appropriate tissue thickness to deliver superior performance in variable tissue.^{3,†}

Poor staple formation is a risk in dense and variable tissues, especially in the colon and rectum following radiation therapy.

Reloads with Tri-Staple™ technology provide consistent performance over a broader range of tissue thicknesses due to its graduated compression profile.^{3,†}

†Pre-clinical results may not correlate with clinical performance in humans.




SPECIFICATIONS




Vascular/medium tissue reload

Staple heights




Medium/thick tissue reload

Staple heights



Extra-thick tissue reload

Staple heights



Reload	Outer row	Middle row	Inner row	Cut line	Inner row	Middle row	Outer row	Intended tissue thickness range
Tan	3.0 mm	2.5 mm	2.0 mm		2.0 mm	2.5 mm	3.0 mm	0.88–1.8 mm
Purple	4.0 mm	3.5 mm	3.0 mm		3.0 mm	3.5 mm	4.0 mm	1.55–2.25 mm
Black	5.0 mm	4.5 mm	4.0 mm		4.0 mm	4.5 mm	5.0 mm	2.25–3.0 mm






SPECIFICATIONS



Grey reload	0.75–1.0 mm (extra-thin to vascular tissue)
Tan reload	0.88–1.8 mm (vascular to medium tissue)
Purple reload	1.5–2.25 mm (medium to thick tissue)
Black reload	2.25–3.0 mm (extra-thick tissue)

Reload Selection Chart



Extra-thin/Vascular reloads contain three rows of 2.0 mm staples (nonvaried height) and contains all of the other features and benefits of an Endo GIA™ reload with Tri-Staple™ technology.

† Preclinical results may not correlate with clinical performance in humans.





SPECIFICATIONS



Vascular/ medium tan reloads	Open	2.0 mm (.79")	2.5 mm (.098")	3.0 mm (.118")
	Closed	0.75 mm (.030")	1.0 mm (.039")	1.25 mm (.049")
	Backspan	3.0 mm (.118")	3.0 mm (.118")	3.0 mm (.118")
	Wire diameter	0.22 mm (.0085")	0.22 mm (.0085")	0.22 mm (.0085")



Medium/thick purple reloads	Open	3.0 mm (.118")	3.5 mm (.138")	4.0 mm (.157")
	Closed	1.25 mm (.049")	1.5 mm (.059")	1.75 mm (.069")
	Backspan	3.0 mm (.118")	3.0 mm (.118")	3.0 mm (.118")
	Wire diameter	0.22 mm (.0085")	0.22 mm (.0085")	0.24 mm (.0094")



Extra-thick black reloads	Open	4.0 mm (.157")	4.5 mm (.177")	5.0 mm (.197")
	Closed	1.75 mm (.069")	2.0 mm (.79")	2.25 mm (.089")
	Backspan	3.0 mm (.118")	3.0 mm (.118")	3.0 mm (.118")
	Wire diameter	0.24 mm (.0094")	0.24 mm (.0095")	0.24 mm (.0096")





REFERENCES

1. Based on internal test report #PCG-007 rev 1, When compared to Ethicon ECHELON FLEX™* green reloads as part of an analysis comparing different stapler designs and their performance and impact on tissues under compression using two-dimensional finite element analysis. September 2 2011.
2. Based on internal test report #2128-002-2, Final analysis of staple line vascularity using MicroCT. When compared to Ethicon ECHELON FLEX™* uniform staple reloads as part of a bench study conducted in murine models, the vascular volume was statistically superior with Tri-Staple™ technology reloads (p = .011). April 27, 2015.
3. Various benchtop testing demonstrating equivalent or superior performance of reloads within a wide range of thicknesses as compared to predicate (510k submission).
4. Based on internal test report #PCG-004, Undercrimp comparisons in increasing pads of foam between Echelon™* and Tri-Staple™ technology (data on file at Medtronic). January 2012.
5. Based on internal test report #PCG-006, Staple Formation Comparison between Medtronic EGIA60AXT and Ethicon ECR60G in an Ex-Vivo Tissue Model (data on file at Medtronic). January 2012.
6. Based on internal test report #PCG-018, 2D FEA of Linear Staplers (data on file at Medtronic). November 2012.
7. Attard JA, Raval MJ, Martin GR, et al. The effects of systemic hypoxia on colon anastomotic healing: an animal model. *Dis Colon Rectum*. 2005 Jul;48(7):1460–1470.



TECHNOLOGY

In the right patient with sufficient space in the pelvis, the radial reload with its innovative profile offers:

- Deeper access in the pelvis in LAR procedures compared to Ethicon CONTOUR™* curved stapler¹
- Better distal margins in rectal resections¹
- Greater maneuverability and visibility in the pelvis¹

Additionally, Tri-Staple™ technology delivers 6 rows of variable height staples and graduated compression that:

- Generates less stress during compression and clamping²
- May allow greater perfusion into the staple line³
- Provides superior performance in variable thicknesses⁴⁻⁸



LESS STRESS

on tissue during compression and clamping²



GREATER PERFUSION

into the staple line³



SUPERIOR PERFORMANCE

in variable tissue thicknesses⁴⁻⁸

> INTRODU

✓ PRODUCT

✓ MISSTE

> ACC

> MOE
DISS

✓ DIVI
STR

> S

> E

>

>

> E

> ANA

> CLO

> OPEN S



REFERENCES



Medtronic



PERFORMANCE

CHALLENGE

SOLUTION

Mobilization and division of the ultra-low/distal rectum is very challenging.

In a cadaver study, four surgeons found the radial reload allowed better visualization and lower access in the pelvis compared with Contour™*¹

Proper distal rectal resection with clear margins are needed to decrease the risk of local recurrence and potentially increase the rate of sphincter salvage in patients with low rectal cancers.¹

In a cadaver study, four surgeons found the radial reload allowed lower access in the pelvis compared with Contour™* in both the coronal and sagittal planes.¹ Distance was measured from the pelvic floor to device.

View study [here](#).

Poor staple formation is a risk in dense and variable tissues, especially in the rectum following radiation therapy.

The extra-thick black radial reload provides the thickest tissue LAR stapler on the market today.⁹


The unique stepped cartridge face and varied height staples approximate tissue to the appropriate tissue thickness to deliver superior performance over a broader range of tissue thickness.⁴⁻⁸



SPECIFICATIONS



Reload Selection Chart



Ordering Information

Description	Cartridge color	Staple size (inner to outer row)
Endo GIA™ radial reload with Tri-Staple™ technology	Purple	3 mm, 3.5 mm, 4 mm
Endo GIA™ radial reload with Tri-Staple™ technology	Black	4 mm, 4.5 mm, 5 mm





SPECIFICATIONS



Competitive analysis: dimensions

Dimensions		
	Endo GIA™ radial reload	Contour™*
Total length	535 mm	406 mm
Functional access length	320 mm	216 mm
Jaw opening	16.3 mm	19 mm
Head width	81 mm	64 mm
Linear staple line (from first staple to end of cut line)	50 mm	38 mm
Full curved staple line	60 mm	48 mm

Best practices:

- The radial reload has a unique design which necessitates a thorough in-service/discussion with each surgeon. This should be scheduled a few days prior to first clinical use, preferably outside of the OR when possible.
- There is a learning curve to using the device — namely becoming familiar with maneuverability options and insertion technique deep in the pelvis. In certain patients with a narrow pelvis, placement may be difficult. In the right patient, the unique geometry of the radial reload may enable the surgeon to get lower than Contour™* in an LAR procedure.¹
- For the first 3–5 cases with the radial reload, be sure your surgeon selects “easier” cases (mid-to upper-rectal tumors in patients with a wider pelvis) that will allow him/her to get comfortable with all the nuances of the device. Once they have gained this experience, they will be fully ready to tackle the more difficult procedures.





REFERENCES

1. Rivadeneira DE, Verdeja JC, Sonoda T. Improved access and visibility during stapling of the ultra-low rectum: a comparative human cadaver study between two curved staplers. *Ann Surg Innov Res*. 2012;6(1):11.
2. Based on internal test report #2128-002-2, Final analysis of staple line vascularity using MicroCT. When compared to Ethicon ECHELON FLEX™* uniform staple reloads as part of a bench study conducted in murine models, the vascular volume was statistically superior with Tri-Staple™ technology reloads (p = .011). April 27, 2015.
3. Based on internal test report #PCG-007 rev 1, When compared to Ethicon ECHELON FLEX™* green reloads as part of an analysis comparing different stapler designs and their performance and impact on tissues under compression using two-dimensional finite element analysis. September 2, 2011.
4. Based on internal test report #PCG-004, Undercrimp comparisons in increasing pads of foam between Echelon™* and Tri-Staple™ technology (data on file at Medtronic). January 2012.
5. Based on internal test report #PCG-006, Staple Formation Comparison between Medtronic EGIA60AXT and Ethicon ECR60G in an Ex-Vivo Tissue Model (data on file at Medtronic). January 2012.
6. Based on internal test report #PCG-018, 2D FEA of Linear Staplers (data on file at Medtronic). November 2012.
7. Based on internal test report #PCG-019, Comparative Test of Medtronic Endo GIA™ Black Reloads with Tri-Staple™ Technology and Ethicon ECHELON FLEX™* Black Reloads (data on file at Medtronic). June 2014.
8. Based on internal test report #PCG-007 rev 1, When compared to Ethicon ECHELON FLEX™* green reloads as part of an analysis comparing different stapler designs and their performance and impact on tissues under compression using two-dimensional finite element analysis. September 2, 2011.
9. Based on internal report # PCG-019, Comparative test of Endo GIA™ black reloads with Tri-Staple™ technology and Ethicon ECHELON FLEX™* black reloads. June 2014.

Endo GIA™ Ultra Universal Stapler

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION
DISSECTION

✓ DIVISION OF
STRUCTURE

> Signia™
stapling system

> Endo GIA™
Tri-Staple™

> Endo GIA™
reload with
technology

> Endo GIA™
Ultra Universal

> ANASTOMOSES

> CLOSURE

> OPEN STEPS

TECHNOLOGY

The Endo GIA™ Ultra Universal stapler provides features such as:

- Three shaft lengths in short, standard, and XL
- Enhanced articulation with 10 points up to 45 degrees
- Forward/backward compatibility with all Endo GIA™ reloads



**ONE-
HANDED
GRASPING
MECHANISM**

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



Medtronic

PERFORMANCE



CHALLENGE

Hand fatigue in manual stapling is challenging, especially when firing in awkward angles on thick tissue during colorectal surgery.

Laparoscopic colorectal surgery is challenging — especially when trying to optimally position the stapler low in the pelvis.

It can be difficult to adequately manage multiple instruments in a laparoscopic procedure.

SOLUTION

Ergonomic design

The Endo GIA™ Ultra Universal stapler has a lightweight handle which fits comfortably in all hand sizes and fires smoothly, even in thick tissue.

Precise articulation

The Endo GIA™ Ultra Universal stapler has up to 45 degree, controlled articulation, which allows for improved access to challenging anatomy.

One-handed grasping mechanism

The Endo GIA™ Ultra Universal stapler has jaws that can be opened and closed with one hand for precise grasping and manipulation of tissue.

Endo GIA™ Ultra Universal Stapler




> INTRODUCTION

∨ PRODUCTS

∨ M



SPECIFICATIONS

Product offering	Reorder code	Shaft length	Description
	EGIAUSHORT	6 cm	Endo GIA™ Ultra Universal short stapler
	EGIAUSTND	16 cm	Endo GIA™ Ultra Universal stapler
	EGIAUXL	26 cm	Endo GIA™ Ultra Universal XL stapler

> CLOSURE

> OPEN STEPS

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



EEA™ Circular Stapler with Tri-Staple™

TECHNOLOGY

Clinical confidence. Because it is proven technology.

The EEA™ circular stapler with Tri-Staple™ technology has advantages over two-row circular staplers because it:

- Puts less stress on tissue during compression and clamping^{1,†}
- Allows greater perfusion into the staple line^{2,†}
- Provides consistent staple performance over a broad range of tissue thicknesses compared to two row circular staplers³⁻⁷

†Preclinical results may not correlate with clinical performance in humans.



LESS STRESS

on tissue during compression and clamping^{1,†}



GREATER PERFUSSION

into the staple line^{2,†}



SUPERIOR PERFORMANCE

in variable tissue thicknesses³⁻⁷

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION
DISSECTION

> DIVISION OF
STRUCTURE

✓ ANASTOMOSES

> EEA™ circular
stapler with Tri-Staple™
technology

> EEA™ staple
Series™ technology

> CLOSURE

> OPEN STEPS



REFERENCES



Medtronic



PERFORMANCE

CHALLENGE

SOLUTION

Ischemia and reduced perfusion may contribute to anastomotic leaks.⁸

EEA™ circular stapler with Tri-Staple™ technology may allow greater perfusion into the staple line compared to two-row circular staplers.^{9,10,†}

Maintaining hemostasis while minimizing risk of ischemia during colorectal surgery is critical to reducing intraoperative blood loss.

EEA™ circular stapler with Tri-Staple™ technology generates less stress on tissue during compression and clamping compared to two-row circular staplers.^{11,†}

Poor staple formation is a risk in dense and variable tissues, especially in the colon and rectum following radiation therapy.

EEA™ circular stapler with Tri-Staple™ technology provides 30 percent additional security to the staple line during the critical healing period compared to two-row circular staplers.¹²⁻¹⁴

†Preclinical results may not correlate with clinical performance in humans.

EEA™ Circular Stapler with Tri-Staple™ Technology

> INTRODUCTION

✓ PRODUCTS

SPECIFICATIONS

Reorder code	Product description	Color	Staple size (inner to outer row)
TRIEEA28MT	EEA™ circular stapler with Tri-Staple™ technology 28 mm medium/thick	Purple	3.0 mm, 3.5 mm, 4.0 mm
TRIEEA31MT	EEA™ circular stapler with Tri-Staple™ technology 31 mm medium/thick	Purple	3.0 mm, 3.5 mm, 4.0 mm
TRIEEAXL33MT	EEA™ circular stapler XL length with Tri-Staple™ technology 33 mm medium/thick	Purple	3.0 mm, 3.5 mm, 4.0 mm

> OPEN STEPS

The proven performance of Tri-Staple™
technology on the EEA™ circular stapler.

TECHNOLOGY PERFORMANCE SPECIFICATIONS



REFERENCES >

Medtronic



REFERENCES

1. Based on internal test report #2128-002-2, Final analysis of staple line vascularity using MicroCT. When compared to Ethicon ECHELON FLEX™* uniform staple reloads as part of a bench study conducted in murine models, the vascular volume was statistically superior with Tri-Staple™ technology reloads ($p = .011$). April 27, 2015.
2. Based on internal test report #PCG-007 rev 1, When compared to ECHELON FLEX™* green reloads as part of an analysis comparing different stapler designs and their performance and impact on tissues under compression using two-dimensional finite element analysis. September 2, 2011.
3. Based on internal test report #PCG-004, Undercrimp comparisons in increasing pads of foam between Echelon™* and Tri-Staple™ technology (data on file at Medtronic). January 2012.
4. Based on internal test report #PCG-006, Staple Formation Comparison between Medtronic EGIA60AXT and Ethicon ECR60G in an Ex-Vivo Tissue Model (data on file at Medtronic). January 2012.
5. Based on internal test report #PCG-018, 2D FEA of Linear Staplers (data on file at Medtronic). November 2012.
6. Based on internal test report #PCG-019, Comparative Test of Medtronic Endo GIA™ Black Reloads with Tri-Staple™ Technology and Ethicon ECHELON FLEX™* Black Reloads (data on file at Medtronic). June 2014.
7. Based on internal test report #PCG-007 rev 1, When compared to ECHELON FLEX™* green reloads as part of an analysis comparing different stapler designs and their performance and impact on tissues under compression using two-dimensional finite element analysis. September 2, 2011.
8. Arrard JA, Raval MJ, Martin GR, et al. The effects of systemic hypoxia on colon anastomotic healing: an animal model. *Dis Colon Rectum*. 2005 Jul;48(7):1460–1470.
9. Based on extrapolation of perfusion studies performed for Endo GIA™ with Tri-Staple™ technology: internal test report #2128-002-2, Final analysis of staple line vascularity using MicroCT. April 27, 2015.
10. Based on internal test report #PCG-007, Media absorbency under clamped conditions. August 6, 2012.
11. Based on internal test report #PCG-30, Comparison of circular staplers: tissue compression profiles as determined by 2-D static axisymmetric finite element analysis (FEA).
12. Based on internal test report #2128-194, Comparison of EEA™ circular stapler with Tri-Staple™ technology to EEA™ circular stapler with DST Series™ technology in colo-colonic and gastro-jejunal anastomoses. August 20, 2015.
13. Based on internal test report #RE0036707, Pilot: comparison of EEA™ circular stapler with Tri-Staple™ technology to EEA™ circular stapler with DST Series™ technology in an esophago-gastrostomy using a canine model. February 25, 2015.
14. Based on internal test report #2128-097, Evaluation of early wound healing events in gastrojejunostomies and colonic anastomosis using a three row EEA™ stapler in canines. August 7, 2013.

EEA™ Stapler with DST Series™

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

> ACCESS

> MOBILIZATION
& DISSECTION

> DIVISION OF
STRUCTURE

∨ ANASTOMOSES

> EEA™ circular
stapler with Tri-Stack™
technology

> EEA™ stapler
with DST Series™ technology

> CLOSURE

> OPEN STEPS

TECHNOLOGY

The EEA™ stapler with DST Series™ technology provides:

- Confidence that comes from more than three decades of clinical experience in colorectal surgery
- Proprietary DST that improves staple formation due to its proprietary rectangular wire cross-section that ensures precise staple formation in challenging applications compared to traditional round wire cross-section¹
- Versatile platform to include two lengths (standard and XL) the most important sizes in colon and rectum



REFERENCES



Medtronic



PERFORMANCE

CHALLENGE

Anastomotic leaks can lead to devastating outcomes. Confidence in proven performance and quality is imperative in this critical firing.

Poor staple formation is a risk in dense and variable tissues, especially in the colon and rectum following radiation therapy.

Access to the anastomotic site can be difficult — especially in more proximal rectal lesions.

SOLUTION

The EEA™ stapler with DST Series™ technology provides reliability, consistency and efficiency from three decades of experience and input from colorectal surgeons around the world.¹

The EEA™ stapler with DST Series™ technology incorporates DST with rectangular wire cross-section that bends more reliably in the intended direction to improve staple formation in challenging applications compared to traditional round wire cross-section.^{1,†}

The EEA™ stapler with DST Series™ technology provides availability of XL length with a 17 percent longer shaft length (35 cm working length) than competitive models for improved access.¹




†Preclinical results may not correlate with clinical performance in humans.



SPECIFICATIONS



Stapler	Product description	Color	Staple size
EEA21	EEA™ stapler with DST Series™ technology 21 mm stapler	Aqua	4.8
EEAXL21	EEA™ stapler with DST Series™ technology XL 21 mm stapler	Aqua	4.8
EEA2135	EEA™ stapler with DST Series™ technology 21 mm stapler with 3.5 mm staples	Aqua	3.5
EEAXL2135	EEA™ stapler with DST Series™ technology XL 21 mm stapler with 3.5 mm staples	Aqua	3.5
EEA25	EEA™ stapler with DST Series™ technology 25 mm stapler	White	4.8
EEAXL25	EEA™ stapler with DST Series™ technology XL 25 mm stapler	White	4.8
EEA2535	EEA™ stapler with DST Series™ technology 25 mm stapler with 3.5 mm staples	White	3.5
EEAXL2535	EEA™ stapler with DST Series™ technology XL 25 mm stapler with 3.5 mm staples	White	3.5
EEA28	EEA™ stapler with DST Series™ technology 28 mm stapler	Blue	4.8
EEAXL28	EEA™ stapler with DST Series™ technology XL 28 mm stapler	Blue	4.8
EEA2835	EEA™ stapler with DST Series™ technology 28 mm stapler with 3.5 mm staples	Blue	3.5
EEAXL2835	EEA™ stapler with DST Series™ technology XL 28 mm stapler with 3.5 mm staples	Blue	3.5
EEA31	EEA™ stapler with DST Series™ technology 31 mm stapler	Green	4.8
EEAXL31	EEA™ stapler with DST Series™ technology XL 31 mm stapler	Green	4.8
EEA33	EEA™ stapler with DST Series™ technology 33 mm stapler	Yellow	4.8
EEAXL33	EEA™ stapler with DST Series™ technology XL 33 mm stapler	Yellow	1.8
EEAORVIL21	EEA™ stapler with Orvil™ device 21 mm (compatible with XL stapler only)	Aqua	N/A
EEAORVIL25	EEA™ stapler with Orvil™ device 25 mm (compatible with XL stapler only)	White	N/A



TRADITIONAL STAPLE

Round wire cross-section is more prone to bending in any direction in challenging applications.¹



DST SERIES™ TECHNOLOGY

Rectangular wire cross-section bends more reliably in the intended direction.¹

EEA™ Stapler with DST Series™ Technology

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

> ACCESS

> MOBILIZATION
DISSECTION

> DIVISION OF
STRUCTURE

∨ ANASTOMOSES

> EEA™ circular
with Tri-Staple
technology

> EEA™ stapler
Series™ technology

> CLOSURE

> OPEN STEPS

REFERENCE

1. Rodeheaver G, Internally sponsored study, Directional Stapling Technology (DST): improved reliability in staple formation in canine models (GIA™ Stapler with DST Series™ technology). 2006.

The EEA™ stapler with DST Series™ technology provides reliable and consistent performance backed by years of clinical experience.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Polysorb™ Braided Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION
DISSECTION

> DIVISION OF
STRUCTURE

> ANASTOMOSES

✓ CLOSURE

> Polysorb™
suture

> Biosyn™
suture


> Maxon™
suture

> V-Loc™ knotless
wound closure device

> OPEN STEPS

TECHNOLOGY

Polysorb™ sutures are composed of Lactomer™* copolymer, which is a synthetic polyester composed of glycolide and lactide (derived from glycolic and lactic acids). They are prepared by coating the suture with a mixture of a caprolactone/glycolide copolymer and calcium stearoyl lactylate.



absorbable suture is needed.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Polysorb™ Braided Absorbable Sutures

PERFORMANCE

CHALLENGE

Wound dehiscence, which is estimated to occur in up to 3.5 percent of patients following surgery.

SOLUTION

The advanced extrusion process of the molecule of Lactomer™* 9-1, the exclusive braiding process, and coating system, give the suture increased strength during the critical wound healing period compared to VICRYL™* sutures,¹⁻⁴ and excellent knot security.^{1,5,6}

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS

REFERENCES

Medtronic

Polysorb™ Braided Absorbable Sutures

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

> ACCESS

> MOBILIZATION
DISSECTION

> DIVISION OF CR
STRUCTURES

> ANASTOMOSIS

∨ CLOSURE


> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device

> OPEN STEPS



SPECIFICATIONS

Polysorb™ sutures are indicated for use in soft tissue approximation or ligation

Tensile strength: 2 weeks: 80% USP
 3 weeks: 30% USP

Absorption profile: 56–70 days

suture with applications in all surgical specialities where a strong absorbable suture is needed.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Polysorb™ Braided Absorbable Sutures

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

> ACCESS

> MOBILIZATION
DISSECTION

> DIVISION OF
STRUCTURE

> ANASTOMOSES

∨ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotted
wound closure

> OPEN STEPS



REFERENCES

1. Based on internal test report #RE00081904 Polysorb benchmarking. 2013.
2. van Ramshorst GH, Nieuwenhuizen J, Hop WC, et al. Abdominal wound dehiscence in adults: development and validation of a risk model. *World J Surg.* 2010 Jan;34(1):20–27.
3. Ammirati CT, Goldman G. Wound closure materials and instruments. *Dermatology. 3rd Edition.* 2012, Elsevier Ltd.
4. Hsu A, Mustoe T. The principles of wound healing. *Plastic Surgery Secrets Plus. 2nd Edition.* 2010, Mosby Inc.
5. Debus ES, Geiger D, Sailer M, Ederer J, Thiede A. Physical, biological and handling characteristics of surgical suture material: A comparison of four different multifilament absorbable sutures. *Eur Surg Res.* 1997;29:52–61.
6. Faulkner B, Gear A, Hellewell T, et al. Biomechanical performance of a braided absorbable suture. *J Long Term Eff Med Implants.* 1996;6(3&4):169–179.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Biosyn™

Monofilament Absorbable Sutures

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

> ACCESS

> MOBILIZATION
DISSECTION

> DIVISION OF
STRUCTURE

> ANASTOMOSES

∨ CLOSURE

> Polysorb™
suture

> Biosyn™
suture


> Maxon™
suture

> V-Loc™ knotless
wound closure device

> OPEN STEPS

TECHNOLOGY

Biosyn™ monofilament absorbable sutures are prepared from Glycomer™* 631, a synthetic polyester composed of glycolide (60%), dioxanone (14%), and trimethylene carbonate (26%).



strength retention and absorption profiles to meet your surgical and tissue healing needs.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Biosyn™ Monofilament Absorbable Sutures

> INTRODUCTION

∨ PRODUCT

∨ MISCELLANEOUS

> APPLICATIONS

> MONOFILAMENT

> DUAL

> ANCHORS

∨ CLINICAL

>

>

>

>

> LOC Knotless
wound closure device

> OPEN STEPS



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

PERFORMANCE

CHALLENGE

Maintaining tensile strength across the wound until tissue tensile strength is adequate.

SOLUTION

Biosyn™ monofilament absorbable suture strength (average knot pull) meets USP and EP specifications.¹



Biosyn™ Monofilament Absorbable Sutures

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

> ACCESS

> MOBILIZATION &
DISSECTION

> DIVISION OF CR
STRUCTURES

> ANASTOMOSIS

∨ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device

> OPEN STEPS



SPECIFICATIONS

Biosyn™ sutures are indicated for use in general soft tissue approximation or ligation

Tensile strength: 2 weeks: 75% USP
 3 weeks: 40% USP

Absorption profile: 90–110 days

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Biosyn™

Monofilament Absorbable Sutures

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

> ACCESS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

∨ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device

> OPEN STEPS

REFERENCE

1. Based on internal report: SUT10181, SUT10182, SUT10183, SUT10184, SUT10185, SUT13059, SUT13060, and SUT13116. PER report Biosyn™ size 1 to 6-0 Bigfoot DQA qualification October 3, 2013. Clinical report Biosyn™ size 1 to 6-0 Bigfoot DQA qualification. February 14, 2014.

strength retention and absorption profiles to meet your surgical and tissue healing needs.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Maxon™ Monofilament Absorbable Sutures

> INTRODUCTION

∨ PRODUCTS

∨ MIS STEPS

> ACCESS

> MOBILIZATION
DISSECTION

> DIVISION OF
STRUCTURE

> ANASTOMOSES

∨ CLOSURE

> Polysorb™
suture

> Biosyn™
suture


> Maxon™
suture

> V-Loc™ knotless
wound closure device

> OPEN STEPS

TECHNOLOGY

Maxon™ and Maxon™ CV monofilament polyglyconate synthetic absorbable sutures (clear or green) are prepared from a copolymer of glycolic acid and trimethylene carbonate.



monofilament absorbable sutures that provide excellent strength over the critical wound healing period.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Maxon™ Monofilament Absorbable Sutures

> INTRODUCTION

▼ PRODUCT

▼ MIS S

> AC

> MC
DIS

> DI
ST

> AN

▼ C

>

>

>

>

wound closure device

> OPEN STEPS

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



PERFORMANCE

CHALLENGE

Short-term absorbable sutures or sutures with insufficient tensile strength can cause post-operative fascial dehiscence and wound complications.

SOLUTION

Maxon™ long-term absorbable sutures tensile strength (average knot pull) meets USP and EP specifications.¹



Medtronic

Maxon™ Monofilament Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

✓ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device

> OPEN STEPS



SPECIFICATIONS

Color:	Green, clear
HPIS code:	755_50_10_10
Latex free:	No
Product style:	Absorbable
Suture length:	24" (60 cm), 30" (75 cm), 36" (90 cm)
UNSPSC:	42312201

the critical wound healing period.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Maxon™ Monofilament Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

✓ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device

> OPEN STEPS

REFERENCE

1. Based on Internal Report Number: RE00007428, RE00007429, SUT09014, SUT09017, SUT09018, SUT100001, SUT10092, SUT10094, SUT13134, SUT13142 & SUT13143. PER & Clinical TMC Supplier change BI to SA DQA Qualification Report Maxon™ Size 3-0 to 7-0 March 30, 2016. PER and Clinical Report Maxon™ Size 1 to 2-0, 3-0 to 6-0 DQA Qualification Report January 28, 2009. PER Report Maxon™ Size 1 to 6-0 DQA Bigfoot Qualification Report April 23, 2014. PER Report Maxon™ D-Tach Size 5-0 DQA Qualification Report May 09, 2014.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS

REFERENCES

Medtronic

V-Loc™ Knotless Wound Closure Device

> INTRODUCTION

✓ PRODUCTS

✓ MIS STEPS

> ACCESS

> MOBILIZATION
DISSECTION

> DIVISION OF
STRUCTURE

> ANASTOMOSES

✓ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture


> V-Loc™ knotless
wound closure device

> OPEN STEPS

TECHNOLOGY

The V-Loc™ device uses a dual-angle cut that creates an anchoring bar while optimizing the strand's integrity.

The V-Loc™ device has a unidirectional barb design that distributes tension across the wound. Meanwhile, circumferential barbs anchor tissue at numerous points, eliminating the need for knots.



closure for patients.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



Medtronic

V-Loc™ Knotless Wound Closure Device

> INTRODUCTION

∨ PRODUCT

∨ MISCELLANEOUS

> ABOUT

> MEDICAL
DISCLAIMER

> DISTRIBUTION
STATEMENTS

> ANATOMY

∨ CLINICAL

>

>

>

>

>

> OPEN STEPS

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



PERFORMANCE



CHALLENGE

Conventional sutures depend on knot integrity for security, and they can potentially cause knot-related complications.

SOLUTION


V-Loc™ knotless wound closure device means no need for knots. There's less risk for the ischemic tissue that can result from tying knots too tight.

wound closure device

closure for patients.



Medtronic



SPECIFICATIONS



	V-Loc™ 90 absorbable device	V-Loc™ 180 absorbable device	V-Loc™ PBT nonabsorbable device
Tensile strength	7 days, 90%; 14 days, 75%	7 days, 80%; 14 days, 75%; 21 days, 65%	Permanent
Absorption profile	90–110 days	180 days	Permanent
Procedural applications	Soft tissue approximation where support is required consistent with the absorption profile	Soft tissue approximation where support is required consistent with the absorption profile	Soft tissue approximation
Color	Undyed, violet	Clear, green	Blue
Composition	Glycolide, dioxanone, and trimethylene carbonate	Copolymer of glycolic acid and trimethylene carbonate	Polybutester
Indications	V-Loc™ 90 device and V-Loc™ 180 absorbable wound closure devices are indicated for soft tissue approximation where use of an absorbable suture is appropriate	V-Loc™ 90 device and V-Loc™ 180 absorbable wound closure devices are indicated for soft tissue approximation where use of an absorbable suture is appropriate	V-Loc™ PBT nonabsorbable wound closure devices are indicated for soft tissue approximation

LigaSure™ Maryland Jaw Device

TECHNOLOGY

Reliable vessel-sealing technology:


- Nonstick nano-coated jaws
- Enhanced blunt dissection compared to straight jaws^{2,†}
- Improved tip visualization compared to straight jaws^{3,‡}

The LigaSure™ Maryland jaw device (LF1944) reduces sticking by more than 38 times compared to the Ethicon ENSEAL™* G2 device.^{4,§}

†23 of 32 surgeons surveyed agreed.

‡31 of 33 surgeons surveyed agreed.

§Tissue sticking to device jaws instances measure over 110 seals per device (ForceTriad™ energy platform).



> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEP

✓ MOBILIZ
DISSEC

> LigaS
jaw de

> LigaS
open

> Sonic
ultras

> DIVISIO
STRUCT

> ANASTO

> CLOSUR



REFERENCES



Medtronic

LigaSure™ Maryland Jaw Device



PERFORMANCE

CHALLENGE

Appropriate skeletonization of IMA is needed to allow confidence when sealing.

Mobilization of the colon requires dissection and vessel sealing which could create the need for frequent instrument exchanges.

SOLUTION

The curved jaws of the LigaSure™ Maryland jaw device allow for easy skeletonization of vessels.^{3,†}

The LigaSure™ Maryland jaw device is a one-step vessel sealer that provides the reliability, control, and consistency² of LigaSure™ technology.

†30 out of 33 surgeons surveyed agreed.



REFERENCES



Medtronic

LigaSure™ Maryland Jaw Device

> INTRODUCTION

∨ PRODUCTS

> M

∨ C



SPECIFICATIONS



	LigaSure™ Maryland jaw device	Ethicon Harmonic ACE+™*
Seal plate length	20.3 mm	13.5 mm
Cut length	18.5 mm	13.5 mm
Jaw aperture	12.7 mm	11.5 mm
Vessel size	7 mm	5 mm
Curved jaw	Yes	Yes
Degrees of curve	22 degrees	17.4 degrees
Shaft lengths	23 cm, 37 cm, 44 cm	14 cm, 23 cm, 36 cm, 45 cm
Ability to cold cut	Yes	No

†Ethicon Harmonic ACE+™* instructions for use.

▪ Cold scissors

TECHNOLOGY PERFORMANCE SPECIFICATIONS



REFERENCES



Medtronic

LigaSure™

Maryland Jaw Device

> INTRODUCTION

∨ PRODUCTS

> MIS STEPS

∨ OPEN STEPS

∨ MOBILIZATION &
DISSECTION

> LigaSure™ Maryland
jaw device

> LigaSure Impact™
open instrument

> Sonicision™ cordless
ultrasonic dissector

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE



REFERENCES

1. Based on LigaSure™ Maryland Jaw Sealer/Divider, Nano-coated [instructions for use]. Boulder, CO: Medtronic; 2016.
2. Based on internal test report #R0035742, Maryland validation labs: Independent surgeon feedback collected during porcine labs in Houston and Los Angeles. April 16 to 18 and April 30 to May 3, 2013.
3. Based on internal test report #RE00071598. Maryland Validations Labs Marketing report. Houston and Los Angeles: independent surgeon feedback collected during porcine labs. April 16–18 and April 30–May 3, 2013.
4. Based on internal test report #RE00073194, Tissue sticking comparison of the Ethicon G2™*. Voyant™* 5 mm Fusion, LigaSure™ LF1737, and LigaSure™ LF1937 devices conducted on porcine tissue using the Force Triad™ energy platform. January 18, 2017.
5. Based on internal test report #RE00140529 rev. A, LigaSure™ Maryland device, nano-coated (LF19XX) tissue testing (memo). March 5, 2018.

TECHNOLOGY PERFORMANCE SPECIFICATIONS



REFERENCES >

Medtronic

LigaSure Impact™

Curved Large Jaw Sealer/ Divider Nano-coated

TECHNOLOGY


Innovative nonstick, nano-coating on large device jaws:

- Results in less cleaning during the procedure than the legacy device (LF4318)^{2,†}
- Makes cleaning more efficient during the procedure than the legacy device (LF4318)^{2,†}

The LigaSure Impact™ device (LF4418) has been shown to have 39 percent fewer less sticking instances than the Ethicon ENSEAL™* G2 Super Jaw.^{3,‡}

†Cleaning effectiveness assessed using porcine uterine tissue, a wet gauze cleaning fixture, and optical imaging analysis after each of 2 cleaning cycles.

‡Tissue sticking to device jaws instances measured over 110 seals per device (ForceTriad™ energy platform).



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS

REFERENCES

Medtronic

LigaSure Impact™

Curved Large Jaw Sealer/

PERFORMANCE

CHALLENGE

Eschar buildup can:

- Cause jaws to stick
- Require multiple cleanings that disrupt procedural flow

Efficient transection of the mesentery requires the ability to take large bites of tissue.

SOLUTION

The LigaSure Impact™ device reduces eschar buildup by 50 percent compared to the legacy device (LF4318), bringing greater efficiency to the operating room (OR).²

The LigaSure Impact™ device (LF4418) enables fast tissue transection.¹

LigaSure Impact™

Curved Large Jaw Sealer/ Divider, Nano-coated

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

✓ MOBILIZATION
DISSECTION

> LigaSure™
jaw device


> LigaSure Impact™
open instr

> Sonicision™
ultrasonic

> DIVISION OF
STRUCTURE

> ANASTOMO

> CLOSURE



SPECIFICATIONS

LigaSure™ technology can seal:

- Vessels up to and including 7 mm
- Lymphatics
- Tissue bundles

Shaft length: 18 cm

Seal length: 36 mm

Cut length: 34 mm

Shaft rotation: 180 degrees

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

LigaSure Impact™

Curved Large Jaw Sealer/ Divider, Nano-coated

> INTRODUCTION

∨ PRODUCTS

> MIS STEPS

∨ OPEN STEPS

∨ MOBILIZATION &
DISSECTION

> LigaSure™ Maryland
jaw device

> LigaSure Impact™
open instrument

> Sonicision™ cordless
ultrasonic dissector

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

> CLOSURE

REFERENCES

1. Based on internal test report #RE00052321, Lig-40 marketing report surgeon evaluation of LigaSure Impact™ device LF4418; independent surgeon feedback collected during Medtronic sponsored porcine labs. June 21–22, 2016.
2. Based on internal benchtop testing #RE00057355, Lig-40 Report LF4418. July 29, 2016.
3. Based on internal benchtop testing #RE00055066, LigaSure™ LF4318, and LigaSure™ LF4418 devices; using porcine uterine tissue, tissue sticking comparison of the Ethicon™* G2 Super Jaw, Olympus ThunderBeat™*. August 2, 2016.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Sonicision™

Cordless Ultrasonic Dissector

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

✓ MOBILIZATION
DISSECTION

> LigaSure™
jaw device

> LigaSure In
open instru

> Sonicision
ultrasonic

> DIVISION OF
STRUCTURE

> ANASTOMO

> CLOSURE

TECHNOLOGY

Offering improved freedom of movement and¹⁻³:

- Innovative cordless technology
- Rapid dissection of tissue and security of hemostasis for vessels up to 5 mm in diameter⁴
- Intuitive energy activation control improves focus on the procedure¹⁻³



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Sonicision™

Cordless Ultrasonic Dissector



PERFORMANCE

CHALLENGE

The growing popularity of laparoscopic and robotic procedures require additional equipment in the OR.

Capital equipment requires planning and the appropriate budgets.

SOLUTION

Sonicision™ cordless ultrasonic device eliminates a cord and it has been shown that the reduction of cords improves safety in the OR.¹

Sonicision™ cordless ultrasonic device platform is easy to adopt without a large capital outlay.

Sonicision™


Cordless Ultrasonic Dissector

> INTRODUCTION

∨ PRO

> M

∨ C



SPECIFICATIONS

Measure	Sonicision™ system	Ethicon Harmonic ACE™†	Ethicon Harmonic ACE+™‡	Olympus Thunderbeat™§
Shaft length	13 cm	14 cm	23 cm	10 cm
	26 cm	23 cm	36 cm	20 cm
	39 cm	36 cm		35 cm
	48 cm	45 cm		45 cm
Method of activation	Hand	Hand	Hand	Hand
	1 button, 2 stages	2 buttons	2 buttons	2 buttons
Frequency	55.5 kHz	55.5 kHz	55.5 kHz	47.7 kHz
Jaw length	14.5 mm	11.1 mm	10.7 mm	16.5 mm
Total weight of device assembled (with cords and transducer)	394.25 g	337.47 g	341.86 g	386.86 g

†Ethicon Harmonic ACE™* instructions for use.
 ‡Ethicon Harmonic ACE+™* instructions for use.
 §Olympus Thunderbeat™* instructions for use.

TECHNOLOGY PERFORMANCE SPECIFICATIONS



REFERENCES >

Sonicision™

Cordless Ultrasonic Dissector

> INTRODUCTION

∨ PRODUCTS

> MIS STEPS

∨ OPEN STEPS

∨ MOBILIZATION & DISSECTION

> LigaSure™ Maryland jaw device

> LigaSure Impact™ open instrument

> Sonicision™ cordless ultrasonic dissector

> DIVISION OF CRITICAL STRUCTURES

> ANASTOMOSIS

> CLOSURE



REFERENCES

1. Based on internal report #R0015050, CCR review for project Milwaukee. Nov. 21, 2013.
2. Based on internal test report #R0014853 Rev A. Usability validation report surgical evaluation of Milwaukee: independent surgeon feedback collected during internally-sponsored cadaver and porcine labs. April 2010 and May 2010.
3. Based on internal report #2-108-10, Validation of Milwaukee USG1, 5 mm cordless laparoscopic ultrasonic shears. April–May, 2010.
4. Tsirline VB, Lau KN, Swan RZ, et al. Evaluation of an innovative, cordless ultrasonic dissector. *Surg Innov.* 2013 Oct;20(5):524–529.
5. Brogmus G, Leone W, Butler L, Hernandez E. Best practices in OR suite layout and equipment choices to reduce slips, trips, and falls. *AORN J.* 2007;86(3):384–394.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

TA™ Stapler with DST Series™ Technology

TECHNOLOGY

DST provides reliability in staple formation during open colorectal surgery:

- Available in 3 staple sizes (2.5 mm, 3.8 mm, and 4.8 mm) to accommodate various tissue thicknesses in the colon and rectum.
- Has a parallel jaw closure with retention pin to ensure consistent B-shaped staples.
- With a strong, low profile anvil, and convenient cutting guide, to optimize access deep in the pelvis.
- Can be reloaded up to 7 times for a total of 8 firings per instrument.



**STRONG,
LOW-PROFILE
ANVIL
WITH
CONVENIENT
CUTTING
GUIDE
OPTIMIZES
ACCESS DEEP
IN THE PELVIS**

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION
DISSECTION

✓ DIVISION OF
STRUCTURE

> TA™ stapler
Series™ tec

> GIA™ stapler
Series™ tec

> SurgiSleeve
wound pro

> ANASTOMO

> CLOSURE



REFERENCES



Medtronic

PERFORMANCE



CHALLENGE

SOLUTION

Poor staple formation is a risk in dense and variable tissues, especially in the colon and rectum following radiation therapy.

Proven DST Series™ staple geometry provides improved superiority in staple strength.¹

Precise stapler positioning can be difficult deep in the pelvis.

The TA™ stapler with DST Series™ technology offers a strong, low profile anvil which provides optimum access without deflection.

Positioning and clamping a stapler on thick tissues can be challenging.


The TA™ stapler has a dual action handle with tactile feedback: first squeeze for closing, second squeeze for firing stapler.



SPECIFICATIONS



TA™ stapler with DST Series™ technology code	Product description	Cartridges available	Staple length before closure	Staple leg length closed	Color	Number of reloads/total number of firings	Number of staples	Staple cross section	Rows of staples
TA30V3S	30 mm	TA30V3L	2.5 mm	1.0 mm	White	7/8	23	.18 x .27 mm	3
TA3035S	30 mm	TA3035L	3.5 mm	1.5 mm	Blue	7/8	11	.19 x .30 mm	2
TA3048S	30 mm	TA3048L	4.8 mm	2.0 mm	Green	7/8	11	.24 x .35 mm	2
TA4535S	45mm	TA4535L	3.5 mm	1.5 mm	Blue	7/8	15	.19 x .30 mm	2
TA4548S	45 mm	TA4548L	4.8 mm	2.0 mm	Green	7/8	15	.24 x .35 mm	2
TA6035S	60 mm	TA6035L	3.5 mm	1.5 mm	Blue	7/8	21	.19 x .30 mm	2
TA6048S	60 mm	TA6048L	4.8 mm	2.0 mm	Green	7/8	21	.24 x .35 mm	2
TA9035S	90 mm	TA9035L	3.5 mm	1.5 mm	Blue	7/8	33	.19 x .30 mm	2
TA9048S	90 mm	TA9048L	4.8 mm	2.0 mm	Green	7/8	33	.24 x .35 mm	2



TRADITIONAL STAPLE

Round wire cross-section is more prone to bending in any direction in challenging applications.¹



DST SERIES™ TECHNOLOGY

Rectangular wire cross-section bends more reliably in the intended direction.¹

TA™ Stapler with DST Series™ Technology

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION
DISSECTION

✓ DIVISION OF
STRUCTURE

> TA™ stapler
Series™ technology

> GIA™ stapler
Series™ technology

> SurgiSleeve
wound protection

> ANASTOMOSIS

> CLOSURE

REFERENCES

1. Rodeheaver G. Internally sponsored study, directional stapling technology (DST): improved reliability in staple formation in canine models (GIA™ Stapler with DST Series™ technology). 2006.

A reloadable linear stapler that places one double staggered row of titanium staples to occlude tissue during open colorectal surgery.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

GIA™ Stapler with DST Series™ Technology

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION
DISSECTION

✓ DIVISION OF
STRUCTURE

> TA™ stapler
Series™ technology

> GIA™ stapler
Series™ technology

> SurgiSleeve
wound protection

> ANASTOMOSES

> CLOSURE

TECHNOLOGY

DST provides reliability in staple formation during open colorectal surgery

Features include:

- Availability in 3 staple sizes (2.5 mm, 3.8 mm, and 4.8 mm) to accommodate various tissue thicknesses in the colon and rectum
- Ability to be reloaded up to 7 times for a total of 8 firings per instrument



TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



PERFORMANCE

CHALLENGE

Poor staple formation is a risk in dense and variable tissues, especially in the colon and rectum following radiation therapy.

Precise stapler positioning can be difficult in tight spaces.

Positioning and clamping a stapler on thick tissues can be challenging.

Cartridges can be "prefired" when handling prior to loading, rendering the cartridge inoperable adding time, and cost to the procedure.

SOLUTION

Proven DST Series™ staple geometry provides improved superiority in staple strength.¹

The GIA™ stapler has an adjustable firing knob giving surgeons the ability to fire from either side of the instrument.

The GIA™ stapler offers a rear or center hinge clamping for one handed positioning.




The GIA™ stapler has a safety lockout which prevents accidental prefiring before clamping.



SPECIFICATIONS



GIA™ stapler with DST Series™ technology	Product description	Cartridges available	Staple length before closure	Approx. staple height closed	Color coded	Number of reloads/ total number of firings	Number of staples	Staple cross section	Rows of staples	Length of staple line
GIA6025S	60 mm	GIA6025L	2.5 mm	1.0 mm	White	7/8	64	.18 x .27 mm	4	64 mm
GIA6038S	60 mm	GIA6038L	3.8 mm	1.5 mm	Blue	7/8	64	.19 x .30 mm	4	64 mm
GIA6048S	60 mm	GIA6048L	4.8 mm	2.0 mm	Green	7/8	64	.24 x .35 mm	4	64 mm
GIA8038S	80 mm	GIA8038L	3.8 mm	1.5 mm	Blue	7/8	84	.19 x .30 mm	4	84 mm
GIA8048S	80 mm	GIA8048L	4.8 mm	2.0 mm	Green	7/8	84	.24 x .35 mm	4	84 mm
GIA10038S	100 mm	GIA10038L	3.8 mm	1.5 mm	Blue	7/8	104	.19 x .30 mm	4	104 mm
GIA10048S	100 mm	GIA10048L	4.8 mm	2.0 mm	Green	7/8	104	.24 x .35 mm	4	104 mm



TRADITIONAL STAPLE

Round wire cross-section is more prone to bending in any direction in challenging applications.¹



DST SERIES™ TECHNOLOGY

Rectangular wire cross-section bends more reliably in the intended direction.¹

GIA™ Stapler with DST Series™ Technology

> INTRODUCTION

∨ PRODUCTS

> MIS STEPS

∨ OPEN STEPS

> MOBILIZATION
DISSECTION

∨ DIVISION OF
STRUCTURE

> TA™ stapler
Series™ technology

> GIA™ stapler
Series™ technology

> SurgiSleeve
wound protection

> ANASTOMOSIS

> CLOSURE

REFERENCES

1. Rodeheaver G. Internally sponsored study, Directional stapling technology (DST): improved reliability in staple formation in canine models (GIA™ Stapler with DST Series™ technology). 2006.

with our proprietary DST

A reloadable linear stapler that places two double-staggered rows of titanium staples and simultaneously cuts and divides tissue during open colorectal surgery.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

TECHNOLOGY

SurgiSleeve™ wound protector provides optimal exposure and wound protection in surgical cases

SurgiSleeve™ wound protectors:

- Are superior in film strength to Alexis™* wound protector/retractors¹
- Offer equivalent or greater exposure versus Alexis™* O with a retraction ring in the large and extra-large sizes as rated by 98 percent of surgeons surveyed^{3,4}
- Are available with large and extra-large retraction rings that are fast to place and remove⁵
- Are faster to roll down compared to Alexis™* O as rated by a majority of surgeons surveyed⁶
- Come in 5 SKUs to cover all your wound protection needs versus 15 SKUs for Alexis™* O wound protector/retractor





PERFORMANCE

CHALLENGE

SOLUTION

Sufficient retraction at incision site.^{5,7-10}

SurgiSleeve™ wound protector has a blue proximal ring that is not only easier to roll down compared to Alexis™* O, but can be rolled down by just one person.†

Adequate exposure to the body.^{5,7-10}

SurgiSleeve™ wound protector has a retraction ring (large and extra-large sizes) providing optimal retraction and visibility compared to Alexis™* O and a stable platform for access to the body.⁶

Potential for wound contamination.^{5,7-10}

SurgiSleeve™ wound protector has three times stronger film material¹ compared to Alexis™* O offering optimal wound protection.

†Survey of 29 General and OB/GYN surgeons of SurgiSleeve™ Wound Protector and Alexis®* O conducted during a product demonstration event October 28–29, 2014 in San Francisco, CA.

SurgiSleeve™ Wound Protector

> INTRODUCTION

∨ PRODUCTS

> MIS STEPS

∨ OPEN STEPS

> MOBILIZATION &
DISSECTION

∨ DIVISION OF CRITICAL
STRUCTURES

> TA™ stapler with DST
Series™ technology


> GIA™ stapler with DST
Series™ technology

> SurgiSleeve™
wound protector

> ANASTOMOSIS

> CLOSURE

Maximized exposure. Wound protection.
Superior strength.¹



SPECIFICATIONS



Available in 5 sizes: extra-small, small, medium, large with retraction ring, and extra-large with retraction ring.

clinically significant reduction in incisional surgical site infections. The study concluded that barrier wound protection of this nature should be considered routine in this type of surgery.²

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



REFERENCES

1. Based on internal report #2151-002, Comparison of mean film material strengths between SurgiSleeve™ wound protector and Alexis™* wound protector, as part of a puncture resistance material evaluation conducted on February 9, 2012. Statistically significant results $p = 0.004$.
2. Reid K, Pockney P, Draganic B, Smith SR. Barrier wound protection decreases surgical site infection in open elective colorectal surgery: A randomized clinical trial. *Dis Colon Rectum*. 2010; 53: 1374–1380.
3. Based on internal report #2151-021, SurgiSleeve™ large wound protector with retraction ring equivalence testing. October 30, 2014, and
4. Based on internal report #2151-023, SurgiSleeve™ extra-large wound protector with retraction ring equivalence testing. November 4, 2014.
5. Voice-of-the-customer survey of 29 general and OB/GYN surgeons of SurgiSleeve™ wound protector and Alexis™* O large and extra-large conducted during a product demonstration event, October 28–29, 2014, in San Francisco, CA.
6. Information in the table was based on (1) Survey of 29 General and OB/GYN surgeons of SurgiSleeve™ wound protector and Alexis™* O large and extra-large conducted during a product demonstration event Oct. 28–29, 2014 in San Francisco, CA; and (2) the animal labs for the large & extra-large ring with five general surgeons and three OB/GYN surgeons on Aug. 5–6 in the North Haven, CT, USA facility.
7. Edwards JP, Ho AL, Tee MC, Dixon E, Ball CG. Wound protectors reduce surgical site infection: a meta-analysis of randomized controlled trials. *Ann Surg*. 2012; 256(1): 53–59.
8. Horiuchi T, Tanishima H, Tamagawa K, et al. A wound protector shields incision sites from bacterial invasion. *Surg Infect (Larchmt)*. 2010;11(6):501–503.
9. Horiuchi T, Tanishima H, Tamagawa K, et al. Randomized, controlled investigation of the anti-infective properties of the Alexis retractor/protector of incision sites. *J Trauma*. 2007;62(1):212–215.
10. Lee P, Waxman K, Taylor B, Yim S. Use of wound-protection system and postoperative wound-infection rates in open appendectomy. *Arch Surg*. 2009;144(9): 872–875.

EEA™ Circular Stapler with Tri-Staple™

TECHNOLOGY

Clinical confidence. Because it is proven technology.

The EEA™ circular stapler with Tri-Staple™ technology has advantages over two-row circular staplers because it:

- Puts less stress on tissue during compression and clamping^{1,†}
- Allows greater perfusion into the staple line^{2,†}
- Provides consistent staple performance over a broad range of tissue thicknesses compared to two row circular staplers³⁻⁷

†Preclinical results may not correlate with clinical performance in humans.



LESS STRESS

on tissue during compression and clamping^{1,†}



GREATER PERFUSION

into the staple line^{2,†}



SUPERIOR PERFORMANCE

in variable tissue thicknesses³⁻⁷

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION
DISSECTION

> DIVISION OF
STRUCTURE

✓ ANASTOMOSES

> EEA™ circular stapler
with Tri-Staple™ technology

> EEA™ stapler
Series™ technology

> CLOSURE



REFERENCES



Medtronic



PERFORMANCE

CHALLENGE

SOLUTION

Ischemia and reduced perfusion may contribute to anastomotic leaks.⁸

EEA™ circular stapler with Tri-Staple™ technology may allow greater perfusion into the staple line compared to two-row circular staplers.^{9,10,†}

Maintaining hemostasis while minimizing risk of ischemia during colorectal surgery is critical to reducing intraoperative blood loss.

EEA™ circular stapler with Tri-Staple™ technology generates less stress on tissue during compression and clamping compared to two-row circular staplers.^{11,†}

Poor staple formation is a risk in dense and variable tissues, especially in the colon and rectum following radiation therapy.

EEA™ circular stapler with Tri-Staple™ technology provides 30 percent additional security to the staple line during the critical healing period compared to two-row circular staplers.¹²⁻¹⁴

†Preclinical results may not correlate with clinical performance in humans.

EEA™ Circular Stapler with Tri-Staple™ Technology

> INTRODUCTION

✓ PRODUCTS



SPECIFICATIONS



Reorder code	Product description	Color	Staple size (inner to outer row)
TRIEEA28MT	EEA™ circular stapler with Tri-Staple™ technology 28 mm medium/thick	Purple	3.0 mm, 3.5 mm, 4.0 mm
TRIEEA31MT	EEA™ circular stapler with Tri-Staple™ technology 31 mm medium/thick	Purple	3.0 mm, 3.5 mm, 4.0 mm
TRIEEAXL33MT	EEA™ circular stapler XL length with Tri-Staple™ technology 33 mm medium/thick	Purple	3.0 mm, 3.5 mm, 4.0 mm

> CLOSURE

The proven performance of Tri-Staple™
technology on the EEA™ circular stapler.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic



REFERENCES

1. Based on internal test report #2128-002-2, Final analysis of staple line vascularity using MicroCT. When compared to Ethicon ECHELON FLEX™* uniform staple reloads as part of a bench study conducted in murine models, the vascular volume was statistically superior with Tri-Staple™ technology reloads ($p = .011$). April 27, 2015.
2. Based on internal test report #PCG-007 rev 1, When compared to ECHELON FLEX™* green reloads as part of an analysis comparing different stapler designs and their performance and impact on tissues under compression using two-dimensional finite element analysis. September 2, 2011.
3. Based on internal test report #PCG-004, Undercrimp comparisons in increasing pads of foam between Echelon™* and Tri-Staple™ technology (data on file at Medtronic). January 2012.
4. Based on internal test report #PCG-006, Staple Formation Comparison between Medtronic EGIA60AXT and Ethicon ECR60G in an Ex-Vivo Tissue Model (data on file at Medtronic). January 2012.
5. Based on internal test report #PCG-018, 2D FEA of Linear Staplers (data on file at Medtronic). November 2012.
6. Based on internal test report #PCG-019, Comparative Test of Medtronic Endo GIA™ Black Reloads with Tri-Staple™ Technology and Ethicon ECHELON FLEX™* Black Reloads (data on file at Medtronic). June 2014.
7. Based on internal test report #PCG-007 rev 1, When compared to ECHELON FLEX™* green reloads as part of an analysis comparing different stapler designs and their performance and impact on tissues under compression using two-dimensional finite element analysis. September 2, 2011.
8. Arrard JA, Raval MJ, Martin GR, et al. The effects of systemic hypoxia on colon anastomotic healing: an animal model. *Dis Colon Rectum*. 2005 Jul;48(7):1460–1470.
9. Based on extrapolation of perfusion studies performed for Endo GIA™ with Tri-Staple™ technology: internal test report #2128-002-2, Final analysis of staple line vascularity using MicroCT. April 27, 2015.
10. Based on internal test report #PCG-007, Media absorbency under clamped conditions. August 6, 2012.
11. Based on internal test report #PCG-30, Comparison of circular staplers: tissue compression profiles as determine by 2-D static axisymmetric finite element analysis (FEA).
12. Based on internal test report #2128-194, Comparison of EEA™ circular stapler with Tri-Staple™ technology to EEA™ circular stapler with DST Series™ technology in colo-colonic and gastro-jejunal anastomoses. August 20, 2015.
13. Based on internal test report #RE0036707, Pilot: comparison of EEA™ circular stapler with Tri-Staple™ technology to EEA™ circular stapler with DST Series™ technology in an esophago-gastrostomy using a canine model. February 25, 2015.
14. Based on internal test report #2128-097, Evaluation of early wound healing events in gastrojejunostomies and colonic anastomosis using a three row EEA™ stapler in canines. August 7, 2013.

EEA™ Stapler with DST Series™

> INTRODUCTION

∨ PRODUCTS

> MIS STEPS

∨ OPEN STEPS

> MOBILIZATION
> DISSECTION

> DIVISION OF
> STRUCTURE

∨ ANASTOMOSES

> EEA™ circular
with Tri-Staple
technology

> EEA™ stapler
with DST Series™ technology

> CLOSURE

TECHNOLOGY

The EEA™ stapler with DST Series™ technology provides:

- Confidence that comes from more than three decades of clinical experience in colorectal surgery
- Proprietary DST that improves staple formation due to its proprietary rectangular wire cross-section that ensures precise staple formation in challenging applications compared to traditional round wire cross-section¹
- Versatile platform to include two lengths (standard and XL) the most important sizes in colon and rectum



REFERENCES



Medtronic



PERFORMANCE

CHALLENGE

Anastomotic leaks can lead to devastating outcomes. Confidence in proven performance and quality is imperative in this critical firing.

Poor staple formation is a risk in dense and variable tissues, especially in the colon and rectum following radiation therapy.

Access to the anastomotic site can be difficult — especially in more proximal rectal lesions.

SOLUTION

The EEA™ stapler with DST Series™ technology provides reliability, consistency and efficiency from three decades of experience and input from colorectal surgeons around the world.¹

The EEA™ stapler with DST Series™ technology incorporates DST with rectangular wire cross-section that bends more reliably in the intended direction to improve staple formation in challenging applications compared to traditional round wire cross-section.^{1,†}



The EEA™ stapler with DST Series™ technology provides availability of XL length with a 17 percent longer shaft length (35 cm working length) than competitive models for improved access.¹

†Preclinical results may not correlate with clinical performance in humans.

SPECIFICATIONS



Stapler	Product description	Color	Staple size
EEA21	EEA™ stapler with DST Series™ technology 21 mm stapler	Aqua	4.8
EEAXL21	EEA™ stapler with DST Series™ technology XL 21 mm stapler	Aqua	4.8
EEA2135	EEA™ stapler with DST Series™ technology 21 mm stapler with 3.5 mm staples	Aqua	3.5
EEAXL2135	EEA™ stapler with DST Series™ technology XL 21 mm stapler with 3.5 mm staples	Aqua	3.5
EEA25	EEA™ stapler with DST Series™ technology 25 mm stapler	White	4.8
EEAXL25	EEA™ stapler with DST Series™ technology XL 25 mm stapler	White	4.8
EEA2535	EEA™ stapler with DST Series™ technology 25 mm stapler with 3.5 mm staples	White	3.5
EEAXL2535	EEA™ stapler with DST Series™ technology XL 25 mm stapler with 3.5 mm staples	White	3.5
EEA28	EEA™ stapler with DST Series™ technology 28 mm stapler	Blue	4.8
EEAXL28	EEA™ stapler with DST Series™ technology XL 28 mm stapler	Blue	4.8
EEA2835	EEA™ stapler with DST Series™ technology 28 mm stapler with 3.5 mm staples	Blue	3.5
EEAXL2835	EEA™ stapler with DST Series™ technology XL 28 mm stapler with 3.5 mm staples	Blue	3.5
EEA31	EEA™ stapler with DST Series™ technology 31 mm stapler	Green	4.8
EEAXL31	EEA™ stapler with DST Series™ technology XL 31 mm stapler	Green	4.8
EEA33	EEA™ stapler with DST Series™ technology 33 mm stapler	Yellow	4.8
EEAXL33	EEA™ stapler with DST Series™ technology XL 33 mm stapler	Yellow	1.8
EEAORVIL21	EEA™ stapler with Orvil™ device 21 mm (compatible with XL stapler only)	Aqua	N/A
EEAORVIL25	EEA™ stapler with Orvil™ device 25 mm (compatible with XL stapler only)	White	N/A



TRADITIONAL STAPLE

Round wire cross-section is more prone to bending in any direction in challenging applications.¹



DST SERIES™ TECHNOLOGY

Rectangular wire cross-section bends more reliably in the intended direction.¹

EEA™ Stapler with DST Series™ Technology

> INTRODUCTION

∨ PRODUCTS

> MIS STEPS

∨ OPEN STEPS

> MOBILIZATION
DISSECTION

> DIVISION OF
STRUCTURE

∨ ANASTOMOSES

> EEA™ circular
with Tri-Staple
technology

> EEA™ stapler
Series™ technology

> CLOSURE

REFERENCE

1. Rodeheaver G, Internally sponsored study, Directional Stapling Technology (DST): improved reliability in staple formation in canine models (GIA™ Stapler with DST Series™ technology). 2006.

The EEA™ stapler with DST Series™ technology provides reliable and consistent performance backed by years of clinical experience.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Polysorb™ Braided Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION
> DISSECTION

> DIVISION OF
> STRUCTURE

> ANASTOMOSES

✓ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device

TECHNOLOGY

Polysorb™ sutures are composed of Lactomer™* copolymer, which is a synthetic polyester composed of glycolide and lactide (derived from glycolic and lactic acids). They are prepared by coating the suture with a mixture of a caprolactone/glycolide copolymer and calcium stearoyl lactylate.



absorbable suture is needed.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Polysorb™ Braided Absorbable Sutures

PERFORMANCE

CHALLENGE

Wound dehiscence, which is estimated to occur in up to 3.5 percent of patients following surgery.

SOLUTION

The advanced extrusion process of the molecule of Lactomer™* 9-1, the exclusive braiding process, and coating system, give the suture increased strength during the critical wound healing period compared to VICRYL™* sutures,¹⁻⁴ and excellent knot security.^{1,5,6}

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS

REFERENCES

Medtronic

Polysorb™ Braided Absorbable Sutures

> INTRODUCTION

∨ PRODUCTS

> MIS STEPS

∨ OPEN STEPS

> MOBILIZATION
DISSECTION

> DIVISION OF CR
STRUCTURES

> ANASTOMOSIS


∨ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device



SPECIFICATIONS

Polysorb™ sutures are indicated for use in soft tissue approximation or ligation

Tensile strength: 2 weeks: 80% USP
 3 weeks: 30% USP

Absorption profile: 56–70 days

suture with applications in all surgical specialities where a strong absorbable suture is needed.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Polysorb™ Braided Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION
> DISSECTION

> DIVISION OF
> STRUCTURE

> ANASTOMOSES

✓ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotted
wound closure



REFERENCES

1. Based on internal test report #RE00081904 Polysorb benchmarking. 2013.
2. van Ramshorst GH, Nieuwenhuizen J, Hop WC, et al. Abdominal wound dehiscence in adults: development and validation of a risk model. *World J Surg.* 2010 Jan;34(1):20–27.
3. Ammirati CT, Goldman G. Wound closure materials and instruments. *Dermatology. 3rd Edition.* 2012, Elsevier Ltd.
4. Hsu A, Mustoe T. The principles of wound healing. *Plastic Surgery Secrets Plus. 2nd Edition.* 2010, Mosby Inc.
5. Debus ES, Geiger D, Sailer M, Ederer J, Thiede A. Physical, biological and handling characteristics of surgical suture material: A comparison of four different multifilament absorbable sutures. *Eur Surg Res.* 1997;29:52–61.
6. Faulkner B, Gear A, Hellewell T, et al. Biomechanical performance of a braided absorbable suture. *J Long Term Eff Med Implants.* 1996;6(3&4):169–179.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Biosyn™

Monofilament Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION
> DISSECTION

> DIVISION OF
> STRUCTURE

> ANASTOMOSES

✓ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device

TECHNOLOGY

Biosyn™ monofilament absorbable sutures are prepared from Glycomer™* 631, a synthetic polyester composed of glycolide (60%), dioxanone (14%), and trimethylene carbonate (26%).



strength retention and absorption profiles to meet your surgical and tissue healing needs.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Biosyn™

Monofilament Absorbable Sutures

> INTRODUCTION

✓ PRODUCT

> MISST

✓ OPEN

> MO
DI

> DI
ST

> AN

✓ cl

> V-Loc™ knotless
wound closure device

PERFORMANCE

CHALLENGE

Maintaining tensile strength across the wound until tissue tensile strength is adequate.

SOLUTION

Biosyn™ monofilament absorbable sutures strength (average knot pull) meets USP and EP specifications.¹

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Biosyn™

Monofilament Absorbable Sutures

> INTRODUCTION

∨ PRODUCTS

> MIS STEPS

∨ OPEN STEPS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

∨ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device



SPECIFICATIONS

Biosyn™ sutures are indicated for use in general soft tissue approximation or ligation

Tensile strength: 2 weeks: 75% USP
 3 weeks: 40% USP

Absorption profile: 90–110 days

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Biosyn™

Monofilament Absorbable Sutures

> INTRODUCTION

∨ PRODUCTS

> MIS STEPS

∨ OPEN STEPS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

∨ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device

REFERENCE

1. Based on internal report: SUT10181, SUT10182, SUT10183, SUT10184, SUT10185, SUT13059, SUT13060, and SUT13116. PER report Biosyn™ size 1 to 6-0 Bigfoot DQA qualification October 3, 2013. Clinical report Biosyn™ size 1 to 6-0 Bigfoot DQA qualification. February 14, 2014.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Maxon™ Monofilament Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION
> DISSECTION

> DIVISION OF
> STRUCTURE

> ANASTOMOSES

✓ CLOSURE

> Polysorb™
suture


> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device

TECHNOLOGY

Maxon™ and Maxon™ CV monofilament polyglyconate synthetic absorbable sutures (clear or green) are prepared from a copolymer of glycolic acid and trimethylene carbonate.



monofilament absorbable sutures that provide excellent strength over the critical wound healing period.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Maxon™ Monofilament Absorbable Sutures

> INTRODUCTION

✓ PRODUCT

> MIS ST

✓ OPEN

> MO
DI

> DI
ST

> AN

✓ cl

>

>

>

>

> Loc-knotless
wound closure device

PERFORMANCE

CHALLENGE

Short-term absorbable sutures or sutures with insufficient tensile strength can cause post-operative fascial dehiscence and wound complications.

SOLUTION

Maxon™ long-term absorbable sutures tensile strength (average knot pull) meets USP and EP specifications.¹

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Maxon™ Monofilament Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS


✓ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device



SPECIFICATIONS

Color:	Green, clear
HPIS code:	755_50_10_10
Latex free:	No
Product style:	Absorbable
Suture length:	24" (60 cm), 30" (75 cm), 36" (90 cm)
UNSPSC:	42312201

the critical wound healing period.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

Maxon™

Monofilament Absorbable Sutures

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION &
DISSECTION

> DIVISION OF CRITICAL
STRUCTURES

> ANASTOMOSIS

✓ CLOSURE

> Polysorb™
suture

> Biosyn™
suture

> Maxon™
suture

> V-Loc™ knotless
wound closure device

REFERENCE

1. Based on Internal Report Number: RE00007428, RE00007429, SUT09014, SUT09017, SUT09018, SUT100001, SUT10092, SUT10094, SUT13134, SUT13142 & SUT13143. PER & Clinical TMC Supplier change BI to SA DQA Qualification Report Maxon™ Size 3-0 to 7-0 March 30, 2016. PER and Clinical Report Maxon™ Size 1 to 2-0, 3-0 to 6-0 DQA Qualification Report January 28, 2009. PER Report Maxon™ Size 1 to 6-0 DQA Bigfoot Qualification Report April 23, 2014. PER Report Maxon™ D-Tach Size 5-0 DQA Qualification Report May 09, 2014.

that provide excellent strength over the critical wound healing period.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



REFERENCES



Medtronic

V-Loc™ Knotless Wound Closure Device

> INTRODUCTION

✓ PRODUCTS

> MIS STEPS

✓ OPEN STEPS

> MOBILIZATION
DISSECTION

> DIVISION OF
STRUCTURE

> ANASTOMOSES

✓ CLOSURE

> Polysorb™
suture

> Biosyn™
suture


> Maxon™
suture

> V-Loc™ knotless
wound closure device

TECHNOLOGY

The V-Loc™ device uses a dual-angle cut that creates an anchoring bar while optimizing the strand's integrity.

The V-Loc™ device has a unidirectional barb design that distributes tension across the wound. Meanwhile, circumferential barbs anchor tissue at numerous points, eliminating the need for knots.



closure for patients.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS



Medtronic

V-Loc™ Knotless Wound Closure Device

> INTRODUCTION

✓ PRODUCT

> MIS ST

✓ OPEN

> MO
DI

> DI
ST

> AN

✓ cl

> V-Loc™ knotless
wound closure device



PERFORMANCE

CHALLENGE

Conventional sutures depend on knot integrity for security, and they can potentially cause knot-related complications.

SOLUTION

V-Loc™ knotless wound closure device means no need for knots. There's less risk for the ischemic tissue that can result from tying knots too tight.

TECHNOLOGY

PERFORMANCE

SPECIFICATIONS






Medtronic



SPECIFICATIONS



	V-Loc™ 90 absorbable device	V-Loc™ 180 absorbable device	V-Loc™ PBT nonabsorbable device
Tensile strength	7 days, 90%; 14 days, 75%	7 days, 80%; 14 days, 75%; 21 days, 65%	Permanent
Absorption profile	90–110 days	180 days	Permanent
Procedural applications	Soft tissue approximation where support is required consistent with the absorption profile	Soft tissue approximation where support is required consistent with the absorption profile	Soft tissue approximation
Color	Undyed, violet	Clear, green	Blue
			
Composition	Glycolide, dioxanone, and trimethylene carbonate	Copolymer of glycolic acid and trimethylene carbonate	Polybutester
Indications	V-Loc™ 90 device and V-Loc™ 180 absorbable wound closure devices are indicated for soft tissue approximation where use of an absorbable suture is appropriate	V-Loc™ 90 device and V-Loc™ 180 absorbable wound closure devices are indicated for soft tissue approximation where use of an absorbable suture is appropriate	V-Loc™ PBT nonabsorbable wound closure devices are indicated for soft tissue approximation