

PROVOX® Vega™ Puncture Set (PVPS)



All-in-one, disposable tracheoesophageal puncture set for primary placement of Provox Vega voice prosthesis

Evidence-based system for both primary and secondary TEPs

Surgeon-Friendly

- Sterile, all-in-one, disposable set ensuring all required instruments are available, and needle is sharp every time
- Minimal preparation and assembly
- Easy dilation and rapid voice prosthesis placement
- 89% immediate insertion with no additional instruments
- 100% successful voice prosthesis placement
- Time reduced by 30% compared to Trocar method
- Extensive onsite clinical and surgical support available

Decreased Time/Cost

- Cost-effective
- Complimentary value analysis assistance provided
- Facility-specific cost comparison of current system vs. Provox Vega Puncture Set

Established Standard of Care

- 265+ US institutions currently using PVPS
- 8,700+ PVPS used in US since introduced in 2014



Provox Vega Puncture Set Disposable Components

Contact your Territory Sales Manager for details!




1 Protect

Pharynx Protector is intended to protect the pharyngeal walls during primary puncture



2 Puncture

Curved, sharp and thin **Puncture Needle** designed for consistently round punctures and minimal tissue trauma



3 Place

The TEP is gradually dilated when pulling the **GuideWire** with the **Dilator** through the opening and the voice prosthesis slides gently in position

Pre-loaded voice prosthesis is designed for a rapid procedure with minimal preparation



Puncture Dilator is designed for gradual dilation to minimize tissue trauma

Retrograde voice prosthesis placement is designed to ensure esophageal flange is placed correctly and allows visual verification of automatic unfolding of tracheal flange

Provox Vega for best voice quality, low speaking effort and ease of cleaning

Atos

Data on File